THUASNE®

GenuPro **ADAPT**

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Conserver cette notice.

en

LIGAMENT AND PATELLA KNEE BRACE WITH SEMI-OPENING AND INTERCHANGEABLE SUPPORTS

Description/Destination

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

Bilateral product.

Available in 7 sizes.

The device is composed of:

- an anatomically-shaped knitting structure including a comfort zone at the back of the knee.
- a semi-opening so it is easier to put on and so it can be adapted to a wider range of thigh morphologies ①.
- an anatomically shaped patella insert ②.
- a pair of removable articulated rigid supports ③.
- a pair of removable flexible stays ④.

This device can be used when resuming sports activities.

Composition

Textile components: polyamide - polyester - elastane - polyurethane - acrylic resin - cotton.

Rigid components: polypropylene - steel - thermoplastic elastomer.

Properties/Mode of action

Elastic fabric with standardized compression ensuring support and proprioception. Anatomically-shaped insert designed to stabilize the patella and contribute to proprioception and pain relief.

Semi-opening on the thigh to make it easier to fit, and so it can be adapted to a wider range of morphologies.

Stabilisation of knee joint ligaments with the rigid hinged side uprights.

Indications

Anterior knee or patella pain.

Mild patellar instability.

Patellar tendinopathy.

Mild to moderate knee osteoarthritis.

Mild joint instability.

Only with articulated rigid supports:

Conservative treatment of minor knee ligament injuries (lateral ligaments).

Contraindications

Do not use the product if the diagnosis has not been confirmed.

Do not use in the event of known allergy to any of the components.

Do not place the product in direct contact with broken skin.

Do not use in the event of major venous thromboembolic history without thromboprophylaxis.

Precautions

Verify the product's integrity before every use.

Do not use the device if it is damaged.

Choose the appropriate size to fit the patient, referring to the size chart.

If the measurements taken do not correspond to the same size, favor the size corresponding to measurement **(b)**.

It is recommended that a healthcare professional supervises the first application.

Strictly comply with your healthcare professional's prescription and recommendations for use.

This product is intended for the treatment of a given condition. Its duration of use is limited to this treatment only.

For hygiene and performance reasons, do not re-use the product for another patient. It is recommended to adequately tighten the device to achieve support without restricting blood circulation.

In the event of discomfort, significant hindrance, pain, variation in limb volume, abnormal sensations or change in colour of the extremities, remove the device and consult a healthcare professional.

In the event of any modification in the product's performance, remove it and consult a healthcare professional.

The use of a medical device by a child should be done under the supervision of an adult or a healthcare professional.

In case of history of venous or lymphatic disorders, consult a healthcare professional.

The lateral uprights have been designed for use with the device only. Do not use the knee brace with other uprights.

Before any sports activity, check the compatibility of the use of this medical device with your healthcare professional.

Do not wear the device in a medical imaging machine.

Do not use the device in case of application of certain products on the skin (creams, ointments, oils, gels, patches...).

Undesirable side-effect

This device can cause skin reactions (redness, itching, burns, blisters...) or wounds of various degrees of severity.

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Application

It is recommended to wear the device directly on the skin, unless contraindicated.

The knee brace is supplied with articulated rigid frames already pre-inserted in the slots provided on each side.

Depending on use, leave them in place or replace them with the flexible stays supplied.

To replace the articulated rigid supports with flexible stays:

On the inside of the knee brace, an opening is located in the upper part of the lateral support covers **1**.

Pull on the tabs to disengage the upper end of the rigid articulated supports from their covers. **2**

Once removed, insert the flexible stays through the cover openings. **3**

Make sure to re-engage the upper end of the covers over the flexible stays.

To replace the flexible stays with articulated rigid supports:

Remove the flexible stays and insert the articulated rigid supports into the covers.

Make sure they are positioned the right way round: the shorter part should be positioned downwards, bending backwards. **4**

Make sure that the curvature of the ends points outwards so it doesn't press against the leg causing discomfort. **5**

To simplify positioning, identification straps have been added to the lower external part of the articulated rigid supports.

Fitting instructions:

To avoid damaging the knitted fabric, make sure that the self-gripping pad is always fastened to the designated gripping surfaces when fitting and removing the brace.

Loosen the strap to allow the device to be semi-opened. **A**

Using the positioning tabs, place the knee brace on the leg, with the foot extended. **B**

Make sure that the patella is positioned in the centre of the patella insert.

Close the strap, adjusting the tightness to the circumference of the thigh. **C**

Check that the brace is pulled completely taut over the leg.

Care

Remove articulated rigid supports before washing, as described in the "Instructions for use" section. Machine washable at 30°C (delicate program). The device has been tested over 20 wash cycles. If possible use a washing net. Do not use detergents, fabric softeners or aggressive products (products containing chlorine, etc.). Do not dry clean. Do not tumble-dry. Do not iron. Squeeze out any excess water. Dry away from any direct heat source (radiator, sun, etc.). If the device is exposed to seawater or chlorinated water, make sure to rinse it in clear water and dry it.

Storage

Store at room temperature, preferably in the original packaging.

Disposal

Dispose of in accordance with local regulations.

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©Thuasne - 2057701 (2025-04)
THUASNE UK Ltd
Thuasne SAS – SIREN/RCS Nanterre 542 091 186