











UniReliever®


fr	Genouillère de décharge, à simple montant articulé	4
en	Offloading single-upright knee brace	5
de	Kniestütze zur Entlastung mit einfacher Gelenkschiene	5
nl	Kniebrace, met enkele scharnierende balein	6
it	Ginocchiera di scarico, con montante articolato semplice	7
es	Rodillera de descarga monoarticular	8
pt	Joelheira de descarga, com uma simples dobradiça articulada	9
da	Aflastende knæbind med en enkelt leddelt stiver	10
fi	Kuormitusta keventävä polvituki yksinkertaisella nivelletyllä pystytuella	10
sv	Knästöd för avlastning, en enkel ledad skena	11
el	Επιγονατίδα αποφόρτισης, με απλή αρθρωτή ενίσχυση	12
cs	Odlehčovací kolenní ortéza s jednou kolenní výztuží	13
pl	Orteza odciążająca z prostą szyną przegubową	14
lv	Celgala ortoze, viens ieliekts statnis	15
lt	Apkrovą mažinantis antkelis su paprasta lankstoma jungtimi	15
et	Koormust vähendav ortoos, üks liigendatud püsttugi	16
sl	Opornica za koleno, enojni pregib	17
sk	Odláhčovacia kolenná ortéza s jednou dlahou	18
hu	Tehermentesítő térdrögzítő, egyoldal ízületes merevítővel	19
bg	Наколенка за разтоварване, едноставен стълб	19
ro	Genunchieră pentru descărcare, cu atelă simplă articulată	20
ru	Наколенник для снятия напряжения с простой шарнирной стойкой	21
hr	Zglobna ortoza za koljeno, s jednostavnom šarkom	22
zh	减压护膝，单侧关节支架	23
ar	داعم ركبة لتخفيف النقل ذو عارضة مفصليّة	24

	6	7
fr	Extension	Flexion
en	Extension	Flexion
de	Streckung	Beugung
nl	Extensie	Flexie
it	Estensione	Flessione
es	Extensión	Flexión
pt	Extensão	Flexão
da	Ekstension	Fleksion
fi	Ekstensio	Fleksio
sv	Extension	Flexion
el	Εκταση	Κάμψη
cs	Extenze	Flexe
pl	Wyprost	Zgięcie
lv	Atliekšana	Saliekšana
lt	Ištiesimas	Sulenkimas
et	Sirutus	Painutus
sl	Iztegovanje	Upogibanje
sk	Extenzia	Flexia
hu	Nyújtás	Hajlítás
bg	Разгъване	Сгъване
ro	Extensión	Flexión
ru	Разгибание	Сгибание
hr	Ekstenzija	Fleksija
zh	伸展	屈曲
ar	المدّ	الثنى

fr	Circonférence de la cuisse (15 cm au-dessus du genou)
en	Thigh circumference (15 cm above the knee)
de	Oberschenkelumfang (15 cm über dem Knie)
nl	Omvang van het dijbeen (15 cm boven de knie)
it	Circonferenza della coscia (15 cm al di sopra del ginocchio)
es	Circunferencia del muslo (15 cm por encima de la rodilla)
pt	Circunferência da coxa (15 cm acima do joelho)
da	Lårets omkreds (15 cm over knæet)
fi	Reiden ympärys (15 cm polven yläpuolelta)
sv	Omkrets runt lår (15 cm ovanför knäet)
el	Περίμετρος του μηρού (15 εκ. πάνω από το γόνατο)
cs	Obvod stehna (15 cm nad kolenem)
pl	Obwód uda (15 cm powyżej kolana)
lv	Augšstilba apkārtmērs (15 cm virs ceļgala)
lt	Šlaunies apimtis (15 cm virš kelio)
et	Reie ümbermõõt (15 cm põlvest kõrgemal)
sl	Obseg stegna (15 cm nad kolenom)
sk	Obvod stehna (15 cm nad kolenom)
hu	Combókörfogat (15 cm a térd felett)
bg	Обиколка на бедрото (15 cm над коляното)
ro	Circumferința coapsei (la 15 cm deasupra genunchiului)
ru	Диаметр окружности бедра (на 15 см выше колена)
hr	Opseg bedra (15 cm iznad koljena)
zh	大腿周长 (膝盖上方15厘米处)
ar	محيط الفخذ (15 سم فوق الركبة)

fr	Décharge	Contrôle du mouvement
en	Off-loading	Motion control
de	Entlastung	Bewegungsführung
nl	Losmaken	Bewegingscontrole
it	Scarico	Controllo del movimento
es	Descarga	Control del movimiento
pt	Alívio do peso	Controlo do movimento
da	Aflastning	Kontrol af bevægelsen
fi	Vapautus	Liikkeen hallinta
sv	Avlastning	Rörelsekontroll
el	Αποφόρτιση	Έλεγχος της κίνησης
cs	Snížení zátěže	Kontrola pohybu
pl	Odciążenie	Kontrola ruchu
lv	Atslogošana	Kustības kontrole
lt	Apkrovos sumažinimas	Judesio kontrolė
et	Koormuse vähendamine	Liikuvuse kontroll
sl	Razbremenitev	Nadzor gibanja
sk	Odlahčenie	Vymedzenie rozsahu pohybu
hu	Tehermentesítés	Mozgásszabályozás
bg	Отбременяване	Контрол на движението
ro	Descărcare	Controlul mișcării
ru	Разгрузка	Контроль движения
hr	Rasterećenje	Kontrola pokreta
zh	减压	运动控制
ar	تخفيف نقل	التحكم في الحركة

fr

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en
OFFLOADING SINGLE-UPRIGHT KNEE BRACE

Description/Destination

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.
Please note that the orthosis can be used as follows:
- Medial arthrosis on the right/lateral arthrosis on the left,
- Medial arthrosis on the left/lateral arthrosis on the right.
One size.

The device is composed of:
- 2 rigid and flexible shells ①,
- 1 lateral upright with TM5+ hinge ② mimicking the knee's natural movement, with protective hinge cap,
- 1 condylar pad,
- 4 inelastic straps (④, ⑤, ⑥, ⑦).
Extension limitations are included in the brace's box (Fig. ①).
Optional flexion limitations ordered separately (Fig. ②).

Composition

Textile components: elastane - polyamide - ethylene vinyl acetate - silicone - polyester - polyurethane.
Rigid components: aluminium - aluminum silicate - acetal - stainless steel - polypropylene - 2-butoxyethanol - dibutyltin oxide - 4-methylpentan-2-one - titanium dioxide - carbon - polyolefin - high density polyethylene - polyamide.

Properties/Mode of action

Adjustable 3-point leverage system with inelastic straps to provide alignment support and off-loading on the affected compartment.
This realignment helps reduce varus or valgus angles of the knee and helps distribute excessive load away from the damaged knee compartment.
Through correction of alignment, the brace will help reduce load, aiding in the reduction of compression and inflammation.
Reducing the load on the affected side would lead to the reduction of pain and cartilage degredation.
The brace maintains proper fit on the leg through:
- inner silicon paddings (shells and straps),
- semi-rigid femoral and tibial shells ① allowing an intimate fit to the patient's unique morphology.
Easy to use thanks to:
- an offloading system with self-dosing adjustment,
- color coded clip buckles.
Extension can be adjusted to 0°, 5°, 10°, 15°, 20° and 30°.
Flexion can be adjusted to 0°, 30°, 45°, 60°, 75°, 90° and 110°.

Indications

Symptomatic unicompartmental femorotibial osteoarthritis (moderate to severe).
Knee off-loading for post-traumatic, postoperative or degenerative conditions.
Joint instability/laxity.
Alternative to osteotomy or leg misalignment surgery.

Contraindications

Do not use the product if the diagnosis has not been confirmed.
Do not use in the event of known allergy to any of the components.
Do not apply the product in direct contact with broken skin or an open wound without an adequate dressing.
Do not use in the event of severe varicose veins (preventing the regular use of an offloading brace).
Do not use in the event of significant genu recurvatum.
Do not use in the event of significant genu varum or genu valgum.
Do not use in the event of major venous thromboembolic history without thromboprophylaxis.

Precautions

Verify the product's integrity before every use.
Do not use the device if it is damaged.
It is recommended that a healthcare professional supervises the first application.
Strictly comply with your healthcare professional's prescription and recommendations for use.
In case of history of venous or lymphatic disorders, consult a healthcare professional.
This product is intended for the treatment of a given condition. Its duration of use is limited to this treatment only.
For hygiene and performance reasons, do not re-use the product for another patient.
It is recommended to adequately tighten the device to achieve support/immobilisation without restricting blood circulation.
In the event of discomfort, significant hindrance, pain, variation in limb volume, abnormal sensations or change in colour of the extremities, remove the device and consult a healthcare professional.
In the event of any modification in the product's performance, remove it and consult a healthcare professional.
Before any sports activity, check the compatibility of the use of this medical device with your healthcare professional.
Do not wear the device in a medical imaging machine.
Do not use the device in case of application of certain products on the skin (creams, ointments, oils, gels, patches...).

Undesirable side-effect

This device can cause skin reactions (redness, itching, burns, blisters...) or wounds of various degrees of severity.
Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Application

It is recommended to wear the device directly on the skin, unless contraindicated.
Please read these instructions with the patient when fitting the brace for the first time and make sure that he/she understands how to position the knee brace.

How to fit the knee brace:

- **Step 1: Evaluate the patient**
Assess the level of pain felt by the patient to measure the function and efficacy of the knee brace after it is fitted.
- **Step 2: Positioning the brace on the leg**
- Unclip and loosen all the straps of the knee brace, and fold the self-fastening straps back for easier application.
- Ensure that the dial is well pulled so that there is no tension in the straps.
- Have the patient sit on the edge of a chair with their knee bent to approximately 30-60° of flexion.
- Place the knee brace on the bare leg (Fig. ③), the hinge on affected compartment side.

Align the hinge at knee center, corresponding with the upper third of your patella (Fig. ③).
- Close the lower calf strap buckle ④, followed by the top thigh strap buckle ⑥ (Fig. ①).
Follow numbers and colors indications on the buckles and the shells.
This should secure the knee brace for positioning of the strap system, ensuring the hinge is positioned correctly.
Position the Cross Buttruss (where the straps intersect) near the edge of patella (Fig. ②).
Use the marking on the Cross Buttruss to ensure proper alignment with the patella (Fig. ②).
- Close front tibia strap buckle ⑤, followed by the front thigh strap buckles ⑥ (Fig. ①).
Match the number and color indications on the buckles to the shells.
- Adjust the back tension of the Cross Straps.
Ensure that the Cross Buttruss is tight against the side of the knee.
Strap chafes can be angled to ensure the straps fit flat against the back side of the leg.

To cut the straps to the desired length, remove the self-fastening hook tab from the end of the strap, cut the strap, and reattach the self-fastening hook tab to the end of the strap (Fig. ②).
Take care to not cut any strap too short.
The comfort pads fixed inside the straps may need to be removed for this operation and repositioned after the strap is cut to avoid any interference when the strap is tightened.
For a better fit to the knee, it is possible to bend the upright to conform to the patient's morphology (Fig. ③):
- Place the knee brace on the edge of a table to stabilize the hinge.
- Bend the upright just above or below the hinge.

• Step 3: Check the fit, suspension and pain level prior to adjusting the corrective force

Ask the patient to stand and walk, taking normal steps and looking straight ahead.
Confirm the brace fits and suspends properly on the leg.
Ask the patient if the pain level with the brace on, still in the neutral position, is the same or less than without the brace.

• Step 4: Adjusting the corrective force (Fig. ⑤)

The dial on the Cross Strap system is set in a neutral position (with dial pulled out).
- To offload the affected side of the knee, push the dial in, and turn clockwise to increase the corrective force to the desired amount based on the level of pain or support needed.
Start with a half turn of the dial and judge whether the correction level still needs to be adjusted after a few steps.
- Turn the dial counter clockwise to lessen the offload or pull to release all the tension.
• **Step 5: Evaluation of pain after adjusting the corrective force**
After adjusting the offloading of the brace:
- Ask the patient to walk again and assess his/her pain level.
- Repeat the operation, increasing or reducing the correction until the proper amount of offloading is achieved for the patient.
- If the patient feels any discomfort, reduce the corrective force.
It is recommended to initially wear the brace for only a few hours a day.
It may take several weeks to feel comfortable with the brace on your leg.

Flexion/extension control:

The brace's default extension setting is 0°.
The extension limitations are attached to a plastic tree in the box (Fig. ④).
• **Adjusting the extension limitation (Fig. ④):**
Extension can be limited to 0°, 5°, 10°, 15°, 20° and 30°.
1. Choose the desired extension limitation on the corresponding plastic piece.
2. Remove the screw located on the hinge side.
3. Flex the hinge slightly and remove the existing extension limitation.
Take care to note the direction this limitation is facing.
4. Insert the desired limitation, numbered end first, with the hook end at the top and facing forward.
Straighten the hinge to full extension to confirm the extension limitation is correctly positioned.
The small hole in each limitation must be aligned and visible through the hole in the hinge, so the screw will thread into the limitation.
5. Reinsert and tighten the screw.

Do a few flexions/extensions to make sure that the limitation is properly locked at the desired angle.

• Adjusting the flexion limitation (Fig. ⑥):

Flexion can be limited to 0°, 30°, 45°, 60°, 75°, 90° and 110°.
1. Choose the desired flexion limitation on the corresponding metallic piece.
Each limitation has the degree etched into the surface.
2. Remove the two screws from the posterior aspect of hinge cover and take out the default spacer.
3. Insert the limitation with the flat end facing up and position it so that both holes in the limitation are visible through the holes in the cap.
4. Insert and tighten the screws through the cap and into both holes in the limitation.

Do a few flexions/extensions to make sure that the limitation is properly locked at the desired angle.

Warning:

The adjustment of the flexion/extension must be decided and performed by the healthcare professional, not by the patient.
NOTE: It is always better to put your brace on a little too high than too low.
Thuasne® cannot be held responsible for undesirable effects or injuries resulting from any unsupervised or inadapted adjustments.

How to take off the knee brace:

To take off the knee brace, pull the dial to release tension and unfasten all straps.

Care

Close the self-fastening tabs before washing. Hand wash. Do not use detergents, fabric softeners or aggressive products (products containing chlorine, etc.). Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.). If the device is exposed to seawater or chlorinated water, make sure to rinse it in clear water and dry it.

• Hinge ②:

The hinge on the brace is pre-lubricated in the factory.
If sand, dirt or water gets inside the hinge, it may require lubricating again.
If you notice the hinge not gliding smoothly, a few drops of a synthetic lubricant can be applied.
Wipe off any excess lubricant before wearing the brace to prevent stains on clothing.

• Straps:

After prolonged use, if the fibers on your strap do not adhere as well to the self-fastening tab, cut the strap shorter so the self-fastening tab adheres to a section of the strap that has fresher fibers.

If this is not possible, you should contact the medical provider who fit your brace.

• Pads:

The brace is lined with padding that provides a comfortable interface between the leg and the shells.

The straps also have pads.

Do not remove these pads from the brace or straps.
Wipe the pads after each use to remove any moisture and let the pads air dry.
You can also clean the pads with a mild anti-bacterial soap and rinse them off with fresh water.
Do not wash the pads in a machine and do not dry in a tumble dryer.

Storage

Store at room temperature, preferably in the original packaging.

Disposal

Dispose of in accordance with local regulations.

COMMERCIAL WARRANTY AGREEMENT AND WARRANTY LIMITATIONS

Thuasne® offers a free, limited commercial warranty to the user, in the territory where the device was purchased, against defects in manufacturing and workmanship for a period of:
- six months for the straps, padding and condylar pad,
- one year for the shells, buckles and the brace's hinge with its cover.

The limited warranty is effective from the date of purchase of the product by the end-user.

The limited commercial warranty does not apply to any defects in manufacturing and workmanship in case of:
- misuse of the product or any damage occurred by a usage outside the normal and intended use of the product as mentioned in the instructions for use,
- damages occurred as part of attempts to modify the product.

Any deterioration or improper trimming of the product during its modification or adjustment process by the healthcare professional when fitting device is expressly excluded from this warranty.

Any claim for this commercial warranty must be sent by the user to the entity where the product was purchased, which will forward this claim to the corresponding Thuasne® entity.

Any warranty claim will first be reviewed by Thuasne® to determine if the conditions of the limited warranty are fulfilled and do not fall into one of the cases of exclusion of the commercial warranty.

To benefit from the warranty, the buyer must mandatorily provide an original and dated proof of purchase of the product.

If the conditions of the limited warranty are fulfilled and the claim is made by the user or its legal representative (parents, guardian...) within the warranty delays indicated above, the buyer will get a new substitution product.

It is expressly agreed that this commercial warranty is in addition to the legal warranties binding the entity which sold the product to the user, in accordance with the applicable local legislation in the country of purchase of the product.

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