THUASNE

PRO MASTER

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Umfang:

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Dear patient

THUASNE has developed the innovative **PROMASTER** bandages as modern orthopaedic aids to help you lead an active and mobile life without discomfort. All therapeutic products in the **PROMASTER** range set standards in respect of shape, styling and functional performance. They are exceedingly comfortable and provide the reassurance you need, even when exercising.

The **Manu PROMASTER** bandage promotes the best possible healing of your injured or painful wrist. Place your trust in this premium bandage featuring the very latest technology.

Manu PROMASTER

Superb mobility. Outstandingly comfortable.

Manu PROMASTER relieves and supports the wrist following injuries, or when patients are suffering from carpal tunnel syndrome or wrist arthrosis. The point-elastic FlexAir structure on the pad channels pressure away from nerves and blood vessels. The soft-knit thumb ensures maximum comfort and mobility.

- Fast and lasting pain reduction
- Supports and stabilises the joint
- Adapts to the wrist without slipping
- Breathable and latex-free knit fabric with 3-D texture guarantees a pleasantly firm fit.

Soft-knit joint

- Maximum comfort at the junction between thumb and index finger

FlexAir pad

- Relieves pressure on the nerves
- Adjusts snugly to the shape of the body
- Alleviates pain

Pull-on aid on the inside arm

- Helps when putting the bandage on

Elasticated compression strap and splint on the inside arm

- Enhanced stabilisation of the wrist

Exercise regime

Wrist

We have developed a short exercise regime for you to promote the healing process. The exercises are designed so that you can perform them anywhere, any time with normal everyday aids. They will take between five and ten minutes, but you can extend or shorten them by altering the number of repetitions. Regular exercise will restore the stability and functional reliability of your wrist more quickly.

EXERCISE AIDS

Chair

2 x 0.5 litre plastic bottles filled with water

NOTES

- Before starting the exercises, make sure you are standing or sitting securely.
- Warm up before exercising; don't stretch cold muscles!
- Stretch gently and carefully just until you begin to feel pain; never go beyond the pain threshold.
- Perform the exercises slowly; focus on what you are doing.
- Continue to breathe calmly and evenly during the exercises.

CAUTION!

None of the exercises should cause pain or exacerbate existing pain. If you have any questions about the exercises themselves, the intensity or frequency, etc., please consult your doctor or physiotherapist.

EXERCISE 1


Warm up – Walking on the spot

Lift your heels one after the other from the floor. Make sure your knees move forward during the lift, not to the side. Try to keep the movements symmetrical on both sides.

Duration:

Walk for 60 seconds – pause for 30 seconds
– walk for 60 seconds

EXERCISE 2

Active stretching – Elbow flexor and extensor


Bend your fingers, wrists and elbows. Turn your forearms inwards while doing so. Hold the position for approx. 10-20 seconds.

Then stretch your fingers, wrists and elbows while turning your forearms outwards. Hold the position for 10-20 seconds.

Duration:

3 repetitions

EXERCISE 3

Stretching the wrist flexor and extensor

Place the palms of your hands together. Bend and stretch the wrists in unison. Keep the movements fluid, don't use force.

Duration:

20 repetitions

EXERCISE 4

Stretching the lateral wrist flexor and extensor


Place the palms of your hands together and bend and stretch the wrists in unison. Keep the movements fluid, don't use force.

Duration:

20 repetitions

EXERCISE 5

Strengthening the arm muscles

Hold your wrist in a neutral position, i.e. hand and forearm form a straight line. Move a light weight (bottles filled with water) back and forth in small rapid movements.

Duration:

Approx. 1 minute

EXERCISE 6

Strengthening – Training the wrist extensor

Lay your forearm on a surface (table) with the inner side facing down. Holding the bottle, let your wrist droop. Stay there for 2-3 seconds before raising your hand by straightening your wrist.

Duration:

10 repetitions on each side,
more if necessary

EXERCISE 7

Strengthening – Training the wrist flexor

Turn your forearm to lay on the outer side without changing the position of your wrist. Holding the bottle, let your wrist droop. Stay there for 2-3 seconds before raising your hand by straightening your wrist.

Duration:

10 repetitions on each side,
more if necessary


MEIN PERSÖNLICHER TRAININGSPLAN MY PERSONAL TRAINING SCHEDULE

Übungen/Exercises:

.....X Wiederholungen/Sets des
Komplettprogramms
Repetitions/sets of the
complete exercise regime

.....X pro Woche/per week:

Meine Trainingstage/my training days

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Recommendations doctor or physiotherapist:

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THUASNE DEUTSCHLAND GmbH