THUASNE®

Genu OA

Genu OA Plus

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Ihr **Thuasne**® Team.

Dear patient

Thuasne® has developed the innovative Genu OA (Plus) knee orthosis as a modern orthopaedic aid to help you lead an active and mobile life.

All therapeutic products made by Thuasne® set standards in respect of shape, styling and functional performance.

They are exceedingly comfortable and provide the reassurance needed in everyday life.

The Genu OA (Plus) orthosis will support the stability and functional reliability of your knee.

Place your trust in this premium orthosis featuring the very latest technology.

Genu OA (Plus)
Innovative,
comfortable,
effective.

Silicone threads in the top knitted edge

- Lessen the threat of the orthosis slipping

Magnetic fasteners at the front

- Make the straps easy to open and close

Cross-over straps

- For dedicated support and optimised pressure distribution

Straps adjusted for anatomical comfort

- Make sure the orthosis sits properly and enhance comfort

Individually adjustable straps on the front and back of the orthosis

- For optimum adjustment and a firm fit

Exercise regime

Knee

We have developed a short exercise regime for you to promote the healing process. The exercises are designed so that you can perform them anywhere, any time with normal everyday aids. They will take between five and ten minutes, but you can extend or shorten them by altering the number of repetitions. Regular exercise will restore the stability and functional reliability of your knee more quickly.

EXERCISE AIDS

Towel

Chair

Cushion

Book


NOTES

- Before starting the exercises, make sure you are standing securely.
- Warm up before exercising: don't stretch cold muscles!
- Stretch gently and carefully just until you begin to feel pain; never go beyond the pain threshold.
- Perform the exercises slowly; focus on what you are doing.
- Continue to breathe calmly and evenly during the exercises.

CAUTION!

None of the exercises should cause pain or exacerbate existing pain. If you have any questions about the exercises themselves, the intensity or frequency, etc., please consult your doctor or physiotherapist.

EXERCISE 1


Warm up – Walking on the spot

Lift your heels one after the other from the floor. Make sure your knees move forward during the lift, not to the side. Try to keep the movements symmetrical on both sides. To increase the intensity, try standing on a cushion to exercise.

Duration:

Walk for 60 seconds – pause for 30 seconds – walk for 60 seconds

EXERCISE 2


Stretching the knee extensor

Bend your supporting leg slightly. Pull your heel towards your bottom using either your hand or a towel. Keep your back straight, avoid a hollow back.

Duration:

Hold for 2 x 15 seconds on each side.

EXERCISE 3


Stretching the rear knee joint

Step back with the affected leg. Both feet are parallel, facing forwards. Move your weight forward and stretch the knee until you can feel the calf muscles stretching in the back leg. Make sure your heel stays on the ground and the back foot doesn't turn outwards.

Duration:

Hold for 2 x 15 seconds on each side.

EXERCISE 4

Strengthening the knee extensor

Sit straight on a chair. Place the moderately heavy book between your feet. Raise the book until your knees are stretched. Hold the position for 5 seconds.

Duration:

5 repetitions

EXERCISE 5

Strengthening the inner thigh muscles

Squeeze a towel or cushion between your knees and press your knees and lower legs together firmly.

Hold for 5-10 seconds before releasing briefly.

Duration:

5 repetitions

EXERCISE 6

Strengthening the knee flexor

Lay on your stomach on a blanket. Place your knees on a thin cushion or rolled-up towel. Rest your head on your folded arms. Bend one knee 90 degrees. Hold the position for about one minute before stretching the leg again. Then change legs.

Duration:

3-5 repetitions on each side

EXERCISE 7

Coordination

Sit upright. Bend your knees at right angles and place both feet on the floor. Raise your heels and lower them again. Then raise your toes and lower them again. Finally, combine both movements: Simultaneously raise the left heel and right toes before lowering them again.

Duration:

5 repetitions

We hope your knee heals quickly and well.


Your **Thuasne®** Team

MEIN PERSÖNLICHER TRAININGSPLAN MY PERSONAL TRAINING SCHEDULE

Übungen/Exercises:

.....X Wiederholungen/Sets des
Komplettprogramms
Repetitions/sets of the
complete exercise regime

.....X pro Woche/per week:

Meine Trainingstage/my training days

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Recommendations doctor or physiotherapist:

NOTIZEN
NOTES

[illegible]

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Sydney Garden 9 · D-30539 Hannover
Tel.: +49 511 20119-0 · Fax: +49 511 20119-177
info@thuasne.de · www.thuasne.com/de-de

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Thuasne Deutschland GmbH