**THUASNE**

**PRO**  
MASTER

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30x

Dear patient

THUASNE has developed the innovative PROMASTER bandages as modern orthopaedic aids to help you lead an active and mobile life without discomfort. All therapeutic products in the PROMASTER range set standards in respect of shape, styling and functional performance. They are exceedingly comfortable and provide the reassurance you need, even when exercising.

The Reversa PROMASTER bandage promotes the best possible healing of your injured or painful lumbar spine. Place your trust in this premium bandage featuring the very latest technology.

## Reversa PROMASTER

### Perfect fit. Reliable support.

Reversa PROMASTER works right where the pain is greatest: The Flex pad guarantees enhanced support and stimulates the lower back muscles – alleviating pain in much the same way as a pressure massage. Thanks to the modified, breathable and softer knit fabric beneath the ribcage, the lumbar bandage fits even more snugly. As an alternative, a silicone pad with nubs can be worn to release tension.

- Alleviates pain and helps speed up the return to mobility
- Massaging pads release tension
- Anatomical bandage design for a snug fit
- Soft edges don't pinch and enhance comfort even further

#### Expansion zones in the ribcage and hip areas

- Bandage does not roll up when seated
- Superbly comfortable

#### Reinforcement rods

- In the back and at the side
- For added support

#### Two grips

- Grips on the inside and outside of the fastening system make it easier to put the bandage on

# Exercise regime

## Lumbar spine

We have developed a short exercise regime for you to promote the healing process. The exercises are designed so that you can perform them anywhere, any time with normal everyday aids. They will take between five and ten minutes, but you can extend or shorten them by altering the number of repetitions. Regular exercise will restore the stability and mobility of your lumbar spine more quickly.

### EXERCISE AIDS

Towel

Chair

Cushion

### NOTES

- Before starting the exercises, make sure you are standing, sitting or lying securely.
- Warm up before exercising: don't stretch cold muscles!
- Stretch gently and carefully just until you begin to feel pain; never go beyond the pain threshold.
- Perform the exercises slowly; focus on what you are doing.
- Continue to breathe calmly and evenly during the exercises.

### CAUTION!

None of the exercises should cause pain or exacerbate existing pain. If you have any questions about the exercises themselves, the intensity or frequency, etc., please consult your doctor or physiotherapist.


### EXERCISE 1

#### Warm up - Walking on the spot

Lift your heels one after the other from the floor. Make sure your knees move forward during the lift, not to the side. Try to keep the movements symmetrical on both sides. To increase the intensity, try standing on a cushion to exercise.

#### Duration:

Walk for 60 seconds – pause for 30 seconds – walk for 60 seconds

### EXERCISE 2

#### Stretching the hip flexor

Keep the leg you want to stretch elongated on the mat. Pull the other leg up towards your chest until you feel the muscles stretching in your groin.

#### Duration:

Hold for 30 seconds on each side.  
2 repetitions

### EXERCISE 3

#### Stretching the hip extensor

Using a towel, pull your heel up towards your bottom until you feel the muscles in your front thigh stretching.

#### Duration:

Hold for 30 seconds on each side.  
2 repetitions

#### EXERCISE 4


#### Stretching the back muscles

Sitting down, let your upper body slowly sink forwards until you feel the muscles stretching in your lower lumbar spine. Breathe into your relaxed back.

**Duration:**

Hold for 30 seconds, 2 repetitions

#### EXERCISE 5

#### Strengthening the abdominal muscles

Place your hands beneath your lower back and press down on them. Continue to breathe calmly.

With legs bent, raise them slightly, one after the other. While doing so, don't lose contact with your hands.

**Make it harder:**

Lift both legs from the mat, one after the other. Keep contact with your hands. Continue breathing.

**Duration:**

30 seconds, 3 repetitions

#### EXERCISE 6


#### Strengthening the back muscles

Press your pelvis into the cushion and then lift your arms and legs from the mat slightly, one after the other. Continue to breathe calmly.

**Duration:**

30 seconds, 3 repetitions

#### EXERCISE 7


#### Strengthening the lateral postural muscles

Raise your hip until your body is in a straight line. Hold for 2 seconds.

**Duration:**

10 repetitions on each side

#### EXERCISE 8

#### Strengthening the upper back muscles and functional postural muscles

Your feet, knees and hips are in a straight line. Bend your knees slightly. Keep your thoracic spine straight; your chin is slightly lowered.

Tense your abdominal muscles lightly and turn your arms outwards; your shoulder girdle slides downwards.

Check your posture in a mirror.

Move your arms gently back and forth together. Make sure your posture doesn't change. Continue breathing.

Move one arm gently back and the other forward, and vice versa. Only go as far and as fast as you can without changing your posture.

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**THUASNE DEUTSCHLAND GmbH**