





THUASNE®


LombaStab Vario

fr	Ceinture de soutien lombaire avec coussin lombaire – Serrage facile et précis**	5
en	Lumbar support belt with lumbar cushion – "Easy, precise tightening"	6
de	LWS-Stabilisierungsorthese mit Lendenkissen – „Einfaches und präzises Justieren“	7
nl	Rugbrace met lumbaalkussen – "Eenvoudige en nauwkeurige aanspanning"	9
it	Cintura di sostegno lombare con cuscino lombare – "Serraggio facile e preciso"	10
es	Faja de sujeción lumbar con almohadilla lumbar de fijación fácil y precisa*	11
pt	Cinta de suporte lombar com almofada lombar – "Aperto fácil e preciso"	12
da	Rygstøttebælte med lændepude – "Nem og præcis stramning"	14
fi	Lannerangan tukivyo ja tyyny – "Helposti ja täsmällisesti kiristettävä"	15
sv	Stödbälte för ryggen med ryggkudde – "Enkel och precis åtdragning"	16
el	Ζώνη οσφυϊκής υποστήριξης με οσφυϊκό μαξιλάρι – «Εύκολη και ακριβής σύσφιξη»*	18
cs	Bederní podpurný pás s bederním polštářem – „Jednoduché a přesné stažení“	19
pl	Podtrzymujący pas lędźwiowy z poduszką lędźwiową – „Łatwe i precyzyjne mocowanie“	20
lv	Muguras jostas daļas atbalsta josta ar polsterējumu – „Viegla un precīza savilkšana“	21
lt	Palaikomas juosmens diržas su juosmens pagalvėle – „Lengvai ir tiksliai suveržiamas“	23
et	Selja- ja nimmepiirkonna tugivöö nimmepadjaga – „Lihtsalt ja täpselt pinguldatav“	24
sl	Križni pas z blazino za ledveni del – „Enostavno in natančno vpetje“	25
sk	Podporný driekový pás s driekovým vankúšom – „Jednoduché a presné dotiahnutie“	26
hu	Ágyéki fűző ágyéki párnával – „Egyszerű és pontos rögzítés“	28
bg	Лумбален колан с опора – „Лесно и точно пристягане“	29
ro	Centură de susținere lombară cu pernă lombară – „Strângere ușoară și precisă“	30
ru	Опорный поясничный корсет с поясничной подушкой – «Легко и просто затянуть»*	32
hr	Lumbalni potporni pojas s jastučićem – „Jednostavno i precizno stezanje“	33
zh	帶有腰墊的腰部穩定護帶-“方便精準收緊”	34
ar	** حزام لدعم أسفل الظهر مزود بوسادة لأسفل الظهر - "شد سهل ودقيق"	37

Instructions d'utilisation - Instructions for use - Gebrauchsanweisung - Gebruiksaanwijzingen - Istruzioni per l'uso - Instrucciones de uso - Instruções de utilização - Anvisningervedrørende brug - Käyttöohjeet - Råd vid användning - Οδηγίες χρήσης - Pokyny k používání - Szczegóły użytkowania - Lietošanas instrukcija - Naudojimo instrukcijos - Kasutusjuhised - Navodila za uporabo - Pokyny na používanie - Használati utasítás - Инструкции за употреба - Instrucțiuni de utilizare - Инструкция по эксплуатации - Upute za uporabu - 使用指南 - إرشادات الاستعمال

	1	60 - 80 cm
	2	75 - 95 cm
	3	90 - 110 cm
	4	105 - 125 cm
	5	120 - 140 cm

fr	Tour de taille
en	Waist circumference
de	Taillenumfang
nl	Tailleomvang
it	Circonferenza vita
es	Contorno de cintura
pt	Perímetro de cintura
da	Taljemål
fi	Vyötäröympärys
sv	Midjemått
el	Περίμετρος μέσης
cs	Obvod pasu
pl	Obwód w talii
lv	Vidukļa apkārtmērs
lt	Juosmens apimtis
et	Vööümbermõõt
sl	Obseg pasu
sk	Obvod pásu
hu	Derékbőség
bg	Обиколка на талията
ro	Circumferința taliei
ru	Обхват талии
hr	Opseg struka
zh	腰围
ar	محيط الخصر

Hauteur : 26 cm

Height: 26 cm

Höhe: 26 cm

Hoogte: 26 cm

Altezza: 26 cm

Altura: 26 cm

Altura: 26 cm

Højde: 26 cm

Korkeus: 26 cm

Höjd: 26 cm

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Výška: 26 cm

Wysokość: 26 cm

Augstums: 26 cm

Aukštis: 26 cm

Kõrgus: 26 cm

Višina: 26 cm

Výška: 26 cm

Magasság: 26 cm

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Conserver cette notice.

en

LUMBAR SUPPORT BELT WITH LUMBAR CUSHION - "EASY, PRECISE TIGHTENING"

Description/Destination

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

Composition

Elastic fabric. ❶

Handles. ❷

Back panel with Coolmax lining. ❸

Lumbar cushion. ❹

Abdominal panels with finger-loops.

Laces.

Textile components: polyamide - polyester - elastane - viscose - polyurethane.

Non-textile components: brushed tempered steel - polypropylene - SEBS - polyamide - polyurethane - polyoxymethylene.

Properties/Mode of action

The device relieves lower back pain by supporting the spine through an increase in intra-abdominal pressure. This reduces the load on the intervertebral discs.

This device uses back pads, implementing a pulley-like system to distribute pressure.

Indications

Non-specific low back pain (acute, subacute and chronic phases).

Lumbar spinal stenosis.

Contraindications

Do not use in the event of known allergy to any of the components.

Do not use on pregnant women.

Do not use in the event of an unstable vertebral fracture.

Do not use in the event of hiatus or abdominal hernia, or of severe liver disease.

Do not use in the event of bone cancer with spinal metastases.

Do not use in the event of circulatory, pulmonary or cardiovascular problems in patients for whom an increase in arterial blood pressure is not recommended.

Precautions

Verify the product's integrity before every use.

Wear the device above thin clothing.

Do not use the device if it is damaged.

Choose the appropriate size to fit the patient, referring to the size chart.

Strictly comply with your healthcare professional's prescription and recommendations for use.

This product is intended for the treatment of a given condition. Its duration of use is limited to this treatment only.

Do not wear the device while sleeping.

In the event of discomfort, significant hindrance, pain, abnormal sensations, remove the device and consult a healthcare professional.

The use of a medical device by a child should be done under the supervision of an adult or a healthcare professional.

Measure the waist circumference of the patient at the mid-distance between the iliac crests and floating ribs.

It is recommended that a healthcare professional supervises the first application.

For hygiene, security and performance reasons, do not re-use the product for another patient.

Do not apply the product in direct contact with skin (risk of pinching).

For safety reasons, do not wear the belt close to machines (risk of becoming caught).

Do not wear the device in a medical imaging machine.

Do not wear the device when driving a vehicle.

Undesirable side-effect

This device can cause skin reactions (redness, itching, irritations...).

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Application

It is recommended that a healthcare professional adjust the stays to the shape of the patient's back before the belt is worn for the first time.

Open the belt and lay it out completely flat so that you can see the handles. ❷

Check that the two tightening handles ❸ are fully retracted.

If necessary, unhook them and separate the two back parts by pulling the ends of the belt, then fasten the handles again.

Stabilisation phase

If recommended by your healthcare professional, position the lumbar cushion ❶ in the belt using the touch fastening system (triangular design pointing down) ❸.

Position the belt around the waist by passing the hands through the hand loops located on the abdominal panels ❹.

The tightening back pads must be positioned outwards (i.e. : not against the skin).

The lacing section must be centered over the spine: the pads should be positioned on either side of the spine.

The bottom of the belt must stop at the top of the crease between the buttocks.

The Thuasne® logo indicates which way up the belt should be put on.

Fasten the belt: right side over left side, without pulling too tight or applying any particular pressure: ❹-❺).

Adjust the tightness using the two handles: pull them forwards then fold them flat against the front of the belt, where they attach like a touch fastener. ❹-❺).

These two handles can be used to adjust the level of compression as required throughout the day ❸.

Mobilisation phase

Remove the lumbar cushion ❶ and use the belt following the same steps as in the immobilisation phase.

Care/Maintenance

Close the self-fastening tabs before washing. Machine washable at 30°C (delicate programme). Do not machine wash the device more than 10 times. Do not dry clean. Do not use detergents, fabric softeners or aggressive products (products containing chlorine, etc.). Do not tumble-dry. Do not iron.

Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.). Dry flat.

Storage

Store at room temperature, preferably in the original packaging.

Disposal

Dispose of in accordance with local regulations.

*Study conducted in-house on a panel of 13 people, July 2019.

This medical device is a regulated health product with CE marking.

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Het Thuasne®-logo duidt de richting van de positie aan. Sluiten van de rugbrace: de rechter zijde bovenop de linker zijde, zonder aan te spannen en zonder een bijzondere druk uit te oefenen (®-©).
Thuasne SAS – SIREN/RCS Nanterre 542 091 186