






THUASNE®

LombaStab Immo

fr	Ceinture de soutien dorso-lombaire avec corset - "Serrage facile et précis"	5
en	Dorso-lumbar support belt with corset - "Easy, precise tightening"	6
de	LWS-Flexionsorthese mit Mobilisierungsfunktion - „Einfaches und präzises Justieren"	7
nl	Dorso-lumbale rugbrace met korsetoptie - "Eenvoudige en nauwkeurige aanspanning"	9
it	Cintura di sostegno dorso-lombare con corsetto - "Serraggio facile e preciso"	10
es	Faja de sujeción dorsolumbar con corsé - "Fijación fácil y precisa"	12
pt	Cinta de suporte dorso-lombar com espartilho - "Aperto fácil e preciso"	13
da	Rygstøttebælte med korset - "Nem og præcis stramning"	15
fi	Selkä- ja lantiotukivyyö korsetilla - "Helposti ja tarkasti kiristettävä"	16
sv	Ryggstöd för ländryggen med korsett - "Enkel och precis åtdragning"	17
el	Ζώνη οσφυϊκής υποστήριξης με κορσέ - «Εύκολη και ακριβής σύσφιξη»	19
cs	Dorso-lumbální podpurný pás s korzetem - „Jednoduché a přesné stažení"	20
pl	Orteza lędźwiowo-krzyżowa z gorsetem - „Łatwe i precyzyjne mocowanie"	22
lv	Muguras krūšu-jostas daļas atbalsta josta ar korseti - „Viegla un precīza savilkšana"	23
lt	Palaikomasis nugaros ir juosmens diržas su korsetu - „Lengvai ir tiksliai suveržiamas"	24
et	Selja- ja nimmepiirkonna tugivöö korsetiga - „Lihtne ja täpne pingutus"	26
sl	Hrbtno-ledvena opornica s steznikom - „Enostavno in natančno vpetje"	27
sk	Podporný chrbtovo-driekový korzet - „Jednoduché a presné dotiahnutie"	28
hu	Háti-ágyéki fűző háttámasztó betéttel - „Egyszerű és pontos rögzítés"	30
bg	Дорзолумбален колан с корсет - „Лесно и точно затягане"	31
ro	Centură de susținere dorso-lombară cu corset - „Strângere ușoară și precisă"	33
ru	Опорный поясничный корсет - «Легко и просто затянуть»	34
hr	Lumbalno-ledni potporni pojas s korzetom - „Jednostavno i precizno stezanje"	36
zh	易于调节的带有胸衣的腰背部稳定护带	37
ar	حزام لدعم الظهر وأسفل الظهر مزود بمشد - "شد سهل"	39

	1	60 - 80 cm
	2	75 - 95 cm
	3	90 - 110 cm
	4	105 - 125 cm
	5	120 - 140 cm

fr	Tour de taille	Hauteur : 35 cm
en	Waist circumference	Height: 35 cm
de	Taillenumfang	Höhe: 35 cm
nl	Tailleomvang	Hoogte: 35 cm
it	Circonferenza vita	Altezza: 35 cm
es	Contorno de cintura	Altura: 35 cm
pt	Perímetro de cintura	Altura: 35 cm
da	Taljemål	Højde: 35 cm
fi	Vyötäröympärys	Korkeus: 35 cm
sv	Midjemått	Höjd: 35 cm
el	Περίμετρος μέσης	Ύψος: 35 εκ.
cs	Obvod pasu	Výška: 35 cm
pl	Obwód w talii	Wysokość: 35 cm
lv	Vidukļa apkārtmērs	Augstums: 35 cm
lt	Juosmens apimtis	Aukštis: 35 cm
et	Vööümbermõõt	Kõrgus: 35 cm
sl	Obseg pasu	Višina: 35 cm
sk	Obvod pása	Výška: 35 cm
hu	Derékbőség	Magasság: 35 cm
bg	Обиколка на талията	Височина: 35 cm
ro	Circumferința taliei	Înălțime: 35 cm
ru	Обхват талии	Высота: 35 см
hr	Opseg struka	Visina: 35 cm
zh	腰围	高度 : 35 cm
ar	محيط الخصر	الارتفاع: 35 سم

fr	Immobilisation	Éducation posturale	Effet proprioceptif
en	Immobilization	Postural education	Proprioceptive effect
de	Ruhigstellung	Erlernen der Körperhaltung	Propriozeptive Wirkung
nl	Immobilisatie	Voorlichting over houding	Proprioceptieve werking
it	Immobilizzazione	Educazione posturale	Effetto propriocettivo
es	Inmovilización	Educación postural	Efecto propioceptivo
pt	Imobilização	Educação postural	Efeito propioceptivo
da	Støtte	Postural genoptræning	Proprioceptisk effekt
fi	Immobilisoi	Asentokuntoutus	Liikeaistiin kohdistuva vaikutus
sv	Immobilisering	Postural träning	Proprioceptiv effekt
el	Ακινητοποίηση	Σωστή στάση σώματος	Αποτέλεσμα ιδιοδεκτικότητας
cs	Znehybnění	Správné držení těla	Proprioceptivní účinek
pl	Unieruchomienie	Korekcja postawy	Propriocepcja
lv	Imobilizācija	Stājas korekcija	Proprioceptīvs efekts
lt	Imobilizavimas	Laikysenos taisymas	Propriocepcinis poveikis
et	Liikumatuks muutmine	Rühitreening	Proprioceptiivne toime
sl	Imobilizacija	Učenje telesne drže	Proprioceptivni učinek
sk	Imobilizácia	Cvičenie posturálnych svalov	Proprioceptívny účinok
hu	Rögzítés	Tartásjavító oktatás	Proprioceptív hatás
bg	Обездвижаване	Приучаване към правилна стойка	Проприоцептивен ефект
ro	Imobilizare	Educație posturală	Efect proprioceptiv
ru	Имобилизация	Исправление осанки	Проприоцептивный эффект
hr	Imobilizacija	Učenje o držanju	Učinak proprioceptivnog vježbanja
zh	固定	姿勢联系	本体感
ar	التثبيت	تدريب وضعي	مفعول استقبال الحس العميق

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Conserver cette notice.

en

DORSO-LUMBAR SUPPORT BELT WITH CORSET - "EASY, PRECISE TIGHTENING"*

Description/Destination

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

Composition

Elastic fabric. ❶

Handles. ❷

Back panel with Coolmax lining. ❸

Abdominal panels with finger-loops. ❹

Thermomouldable corset. ❺

Laces.

Textile components: polyamide - polyester - elastane - viscose - polyurethane.

Non-textile components: brushed tempered steel - polypropylene - SEBS - polyamide - polyurethane - polyoxymethylene - polyethylene.

Properties/Mode of action

The device relieves dorso-lumbar pain by providing support to the lumbar spine.

During the immobilisation phase, firm support is provided thanks to a rigid thermomouldable corset.

During the mobilisation phase, compartmental lumbar unloading is achieved by increasing the intra-abdominal pressure coupled with a 3-pressure-point system.

This device uses back pads, implementing a pulley-like system to distribute pressure.

Indications

Post-operative immobilisation.

MODIC type 1 discopathy.

Non-specific low back pain (acute, subacute and chronic phases).

Spondylolysis.

Spinal stenosis.

Contraindications

Do not use in the event of known allergy to any of the components.

Do not use on pregnant women.

Do not use in the event of an unstable vertebral fracture.

Do not use in the event of hiatus or abdominal hernia, or of severe liver disease.

Do not use in the event of bone cancer with spinal metastases.

Do not use in the event of circulatory, pulmonary or cardiovascular problems in patients for whom an increase in arterial blood pressure is not recommended.

Precautions

Verify the product's integrity before every use.

Wear the device above thin clothing.

Do not use the device if it is damaged.

Choose the appropriate size to fit the patient, referring to the size chart.

Strictly comply with your healthcare professional's prescription and recommendations for use.

This product is intended for the treatment of a given condition. Its duration of use is limited to this treatment only.

Do not wear the device while sleeping.

In the event of discomfort, significant hindrance, pain, abnormal sensations, remove the device and consult a healthcare professional.

The use of a medical device by a child should be done under the supervision of an adult or a healthcare professional.

Measure the waist circumference of the patient at the mid-distance between the iliac crests and floating ribs. It is recommended that a healthcare professional supervises the first application.

For hygiene, security and performance reasons, do not re-use the product for another patient.

Do not apply the product in direct contact with skin (risk of pinching).

For safety reasons, do not wear the belt close to machines (risk of becoming caught).

Do not wear the device in a medical imaging machine.

Do not wear the device when driving a vehicle.

Undesirable side-effect

This device can cause skin reactions (redness, itching, irritations...).

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Application

It is recommended that a healthcare professional adjust the stays to the shape of the patient's back before the belt is worn for the first time.

Open the belt and lay it out completely flat so that you can see the handles. ②

Check that the two tightening handles ② are fully retracted.

If necessary, unhook them and separate the two back parts by pulling the ends of the belt, then fasten the handles again.

Immobilisation phase

Remove the corset cover ④ and fit the latter to the patient's anatomy using a heat source. This must be done by a healthcare professional. This operation must be done by a healthcare professional.

Once cool, return the formed corset to its cover and position this latter in the belt using the self-fastening system (widest part at bottom). ④

Position the belt around the waist, placing your hands in the hand loops located on the front panels. ⑥

The corset must be inside the belt and the tightening back pads ① must be positioned outwards (i.e. : not against the skin).

The lacing section must be centered over the spine: the pads should be positioned on either side of the spine. ⑥

The bottom of the belt must stop at the top of the crease between the buttocks. ⑥

The Thuasne® logo indicates which way up the belt should be put on.

Fasten the belt: right side over left side, without pulling too tight or applying any particular pressure. ⑥ ⑦

Adjust the tightness using the two handles: pull them forwards then fold them flat against the front of the belt, where they attach like a self-fastener. ⑥ ⑦

These two handles can be used to adjust the level of compression as required throughout the day.

Mobilisation phase

Remove the thermomouldable corset that gives the belt its rigidity ①, then use the belt in the same manner as during the immobilisation phase.

Care

Corset :

Clean the corset with a damp cloth.

Belt and cover:

Close the self-fastening tabs before washing.

Machine washable at 30°C (delicate programme).

Do not machine wash the device more than 10 times. Do not dry clean. Do not use detergents, fabric softeners or aggressive products (products containing chlorine, etc.). Do not tumble-dry. Do not iron. Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.). Dry flat.

Storage

Store at room temperature, preferably in the original packaging.

Disposal

Dispose of in accordance with local regulations.

*Study conducted in-house on a panel of 13 people, July 2019.

This medical device is a regulated health product with CE marking.

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©Thuasne - 2054301 (2024-04)
THUASNE UK Ltd
Thuasne SAS – SIREN/RCS Nanterre 542 091 186