

THUASNE®



VENOFLEX

de Gebrauchsanweisung - Medizinische Kompressionsstrümpfe
en Instructions for use - Medical compression stockings or tights
fr Consignes d'utilisation - Bas ou collants de compression médicaux
it Istruzioni per l'uso - Calze e collant di contenzione terapeutica
es Instrucciones de utilización - Medias de compresión terapéuticas
pt Instruções de uso - Meias de compressão medicinais
nl Gebruiksaanwijzing - Medische compressiekousen
da Brugsanvisning - Medicinske kompressionsstrømper
sv Bruksanvisning - Medicinska kompressionsstrumpor
fi Käyttöohje - Lääketieteelliset puristussukat
no Bruksanvisning - Medisinske kompresjonsstrømper
el Οδηγίες χρήσης - Ιατρικές κάλτσες διαβαθμισμένης συμπίεσης
cs Návod k použití - Zdravotní kompresivní punčochy
pl Sposób użytkowania - Medyczne wyroby uciskowe
et Kasutusjuhend - Meditsiinilised kompressioonsukad
lv Lietošanas instrukcija - Medicīniskās kompresijas zeķes
lt Naudojimo instrukcija - Medicininės kompresinės kojinės
sk Návod na použitie - Zdravotné kompresívne pančuchy
sl Navodila za uporabo - Medicinske kompresijske nogavice
hu Használati utasítás - Gyógyászati kompressziós harisnya
ro Instrucțiuni de utilizare - Ciorapi medicinali compresivi
bg Инструкция за употреба - Медицински компресивни чорапи
hr Upute za uporabu - Medicinske kompresijske čarape
ru Инструкция по использованию - Медицинские компрессионные чулки
vi Hướng dẫn sử dụng - Vớ quần áp lực y tế
ko 독일어 - 의료용 압박 스타킹
ja 使用説明書 - 医療用弾性ストッキング
zh 使用說明 - 医用彈力襪
ar تعليمات الاستخدام - جوارب ضغط طبية


1.

2.

3.

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Purpose

Compression stockings are the basic therapy in the treatment of venous disorders and oedemas. In addition to their medical effect, compression stockings support your well-being as a whole. Compression therapy comprises measures of exterior, planar pressure applications in case of venous disorders and lymphatic drainage problems. The application of planar pressure is to prevent the development of oedemas and support venous reflux or lymphatic drainage. Compression therapy generally is a long-term treatment and requires the corresponding products which – if the medical indication is given – the patient him/herself applies upon the doctor's instruction.

Indications

Primary and secondary varicosis; varices during pregnancy; supporting sclerotherapy; after surgical venous interventions; superficial thrombophlebitis, as well as condition after healed up phlebitis; deep vein thrombosis; condition after thrombosis; postthrombotic syndrome; thrombosis prophylaxis for mobile patients; CVI of stages I – III according to Widmer and C1 – C6 according to CEAP; ulcer prevention; venous insufficiency; angiodysplasias; lymphoedemas; oedemas during pregnancy; post-traumatic oedemas; postoperative oedemas; cyclic idiopathic oedemas; lipoedemas; lipolymphoedemas; phlebo-lymphatic oedemas; congestive conditions as a result of immobility (arthrogenous congestion syndrome; pareses and partial pareses of extremities); improvement of venous symptoms; occupational oedemas (walking and standing work; sedentary work). After burns and for scar treatment, we recommend our special products for scar therapy.

Absolute contraindications

Peripheral arterial occlusive disease (PAOD) stage III/IV; decompensated heart failure; septic phlebitis; phlegmasia coerulea dolens; acute erysipelas; acute bacterial, viral or allergic inflammations; untreated or advanced dermatoses /contact dermatitis; untreated open wounds; swelling of extremities of unknown origin; untreated malignant tumours.

Relative contraindications

Pronounced exudative dermatoses; malignant lymphoedema; peripheral arterial occlusive disease (PAOD) stage I/II; peripheral neuropathy associated with diabetes mellitus; sensitivity disorders of extremities; compensated heart failure; skin infections; dermatoses (e.g. blistering dermatoses); intolerance of pressure or of product ingredients; primary chronic polyarthritis. In case of open wounds in the treatment area, the wound area must first be covered with an adequate wound dressing or a bandage before the compression stockings are put on.

Thuasne Deutschland GmbH accepts no liability for failure to observe these contraindications.

Side effects

No side effects are known if properly used. In case of any negative changes during the application, please consult your doctor or the medical retailer immediately. In case of complaints in connection with the product, such as damage of the knitted fabric or lack of fit, please contact directly your medical retailer.

The manufacturer shall not be liable for damages / injuries which resulted due to improper handling or misuse.

Precautions

Follow the instructions given by your doctor, by an orthopaedic technician or medical expert personnel. If you experience severe pain or discomfort, please see your doctor, an orthopaedic technician or medical expert personnel. Store at room temperature, preferably in the original packaging.

To guarantee its effectiveness and for reasons of hygiene, the product may not be passed on to other patients. Any serious incident that occurs in connection with the product must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient are/is domiciled.

Directions for putting the stockings on

Put on your compression stockings in the morning after getting up. Sharp finger nails or jewellery may damage the stockings. It would be best to put on your compression stockings with rubber gloves; you can thus roll up the stocking with the palm of your hand until it fits perfectly without wrinkles. When you use stockings with an adhesive band, please be sure not to get your

hands into the seam of the adhesive band when putting the stockings on or taking them off.

1. Turn over the compression stocking up to the foot part and pull the foot part of the stocking over the foot. Up to the heel, the foot part must fit well on the foot. For stockings without a toe tip, it is recommended to first put the enclosed foot slide over the foot.

2. Now the compression stocking can be carefully rolled up the leg until it is pulled on to its full length. The stocking should not be pulled up forcefully or overstretched by grasping the top edge of the stocking or its adhesive band.

3. To make the stocking fit well, push up the stocking material in „waves“ to the top beginning at the bottom – if necessary, do it in several steps; and finally, correct the fit with the palm of your hand.

The stocking may be easily damaged due to rough and cracked skin on heels. Thus make sure, when putting the stocking on, that your heel gently moves through the tight leg part of the stocking.

Please also ensure that the inside of your footwear is perfect. Defective lining of footwear increases friction and may damage the stocking.

We recommend our fitting aids if there are difficulties putting on the stockings.

Care

In order to maintain their function, wash your stockings every day after their removal. Use a mild detergent for washing by hand or for machine wash (gentle cycle at 40° C) without optical brightener and fabric softener. In a laundry bag, compression stockings are specially protected. In order to reduce the drying period, put the compression stockings into a terry cloth towel, roll it up firmly and press the stockings out (do not wring them out). Do not leave the moist stockings in the wet towel. Dry your stockings in the low cycle of your dryer or simply let them air dry. Never put them in the sun or onto a radiator to dry. Do not expose the stockings to any damaging effects, such as: strong sunlight, intense heat, ointments, oils and greases, gasoline and aggressive detergents.

Material composition

Please find the exact information on the sewn-in fabric label.

Note: The sewn-in fabric label is an important component of the compression stockings. If you remove this labelling from a Thuasne® product, any claim for warranty, repair or replacement shall be waived.

Storage information and duration of use

Store in a dry place and protect from sunlight. For series products, a maximum shelf life of 36 months applies. This is followed by the duration of use of 6 months for the compression clothing. The hourglass pictogram shows the duration of use printed on the carton label. Customised manufacture concerns products for immediate use.

We guarantee the medical efficacy in series sizes and customised manufacture for a maximum wearing period of six months. This is subject to the correct handling of the stockings (for example, their care or putting them on or taking them off). If necessary, your doctor should then prescribe new compression stockings. In case of any significant change of relevant body sizes (e.g. due to a successful therapy, weight change), new measurements and a new supply may also be justified earlier already.

Disposal

You may dispose of your medical compression stockings via the residual or household waste. No special criteria exist for their elimination.

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Thuasne Cervitex Ltd. · The Bell House, 20 Hataas St. · Kfar Saba, 44425
Thuasne Scandinavia AB · Box 95 · 131 07 Nacka
Thuasne UK Ltd