
BIPLAST

**fr** Bande adhésive de contention  
**en** Adhesive compression bandage  
**de** Haftender Kompressionsverband  
**nl** Zelfklevend steunverband  
**it** Fascia adesiva di contenimento  
**es** Venda adhesiva de contención  
**pt** Banda adesiva de contenção  
**fi** Liimasiide  
**sv** Självhäftande stödband  
**el** Αυτοκόλλητη ταινία σύγκράτησης  
**pl** Pas podtrzymujący samoprzylepny  
**lv** Liposa balsta saite  
**sk** Adhezívny fixačný obváz  
**hu** Öntapadós rugalmas fasíl  
**bg** Адезивен бинт за захранка  
**ro** Banda adeziva de compresie  
**ru** Бинт липкий фиксирующий  
**ko** 접착식 고정 봉대  
**ar** بند ملتصق لاصق

**fr** BANDE ADHÉSIVE DE CONTENTION

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Conserver cette notice.

**en ADHESIVE COMPRESSION BANDAGE****Description/Destination**

Short-stretch adhesive compression bandage.  
The device is available in various widths according to needs.  
The device is intended for use on the patient's body part or joint affected by this condition, and only for the treatment of the indications listed.

**Composition**

Coton, polyester, -adhésive -mass (synthetic - elastomer and -tackifying resin), paper.

**Properties/Mode of action**

Adhesive compression bandage.

Aerated elastic fabric.

Multi-directional herringbone coating.

Protective, release backing paper.

The device can be used to secure dressings and splints.

**Indications**

Soft tissue or musculoskeletal system injuries.

Strapping: post-trauma support.

Small bone fractures (finger, foot, etc.).

Congenital or acquired bone deformity of the foot.

Vascular or lymphatic conditions requiring a short-stretch or multicomponent bandaging.

**Contraindications**

Do not use the product if the diagnosis has not been confirmed.

Do not use in the event of known allergy to any of the components.

Do not apply the product in direct contact with broken skin or an open wound without an adequate dressing.

Do not use to keep dressings in place on damaged skin or following an operation.

Do not use in the event of oozing or severe dermatological conditions of the limb.

Do not use in the event of extra-anatomical bypass of the affected limb.

Do not use in the event of severe peripheral neuropathy of the affected limb.

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**Precautions**

Verify the product's integrity before every use.

Do not use the device if it is damaged.

Strictly comply with your healthcare professional's prescription and recommendations for use.

In case of particularly sensitive skin, follow the advice of a healthcare professional.

In the event of any modifications in the product's performance, remove it and consult a healthcare professional.

In the event of discomfort, significant hindrance, pain, variation in limb volume, deterioration of skin condition, infection, abnormal sensations or change in the colour of the extremities, remove the device and seek the advice of a healthcare professional.

Do not use the device in case of application of certain products on the skin (creams, ointments, oils, gels, patches...).

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Pull on both side to tear the full length ③.

Peel off the paper on one side, apply on the limb, then peel off on the other side ④.

**Storage**

Store at room temperature, in the original packaging.

**Disposal**

Dispose of in accordance with local regulation.

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