# THUASNE

## Ligaflex®

- fr** Sangles claviculaires
- en** Clavicular straps
- de** Schlüsselbeingurte
- nl** Sleutelbeenbanden
- it** Immobilizzatore clavicolare
- es** Correas claviculares
- pt** Correias claviculares
- da** Kravebenstropper
- fi** Solisluiden hihnat
- sv** Nyckelbensrem
- el** Ιμάντες κλείδας
- cs** Klavikulární bandáž
- pl** Pasy unieruchamiające obojczyk
- lv** Atslēgas kaula siksnas
- lt** Raktikaulių diržai
- et** Rangluud fikseerivad rihmad
- sl** Trakovi za ključnice
- sk** Klavikulárne popruhy
- hu** Kulcsosontrögzítő
- bg** Ремъци за ключици
- ro** Curele pentru claviculă
- ru** Ремни для ключицы
- hr** Trake za ključne kosti
- zh** 锁骨束带
- ar** أحزمة للترقوة

Instructions d'utilisation - Instructions for use  
 Gebrauchsanleitung - Gebruiksaanwijzingen  
 Istruzioni per l'uso - Instrucciones de uso  
 Instruções de utilização - Anvisningar vedrørende brug - Käyttöohjeet - Råd vid användning  
 Οδηγίες χρήσης - Pokyny k používání - Szczegóły użytkowania - Lietošanas instrukcija - Naudojimo instrukcijos - Kasutusjuhised - Navodila za uporabo  
 Pokyny na používanie - Használati utasítás  
 Инструкции за употреба - Instrucțiuni de utilizare  
 Инструкция по эксплуатации - Upute za uporabu  
 使用指南 - إرشادات الاستعمال

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en  
**CLAVICULAR STRAPS**

**Description/Destination**

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

Figure of 8 straps allowing stabilisation of the clavicle for post-traumatic use.

**Composition**

Cotton - Polyester - Elastane - polyamide - Polyoxymethylene - Ethylene vinyl acetate - Polyurethane.

**Properties/Mode of action**

Two dorsal straps fastened onto an anchoring pad. The simultaneous tightening of the straps helps stabilise the clavicle.

Adjustable and trimmable straps according to patient's shoulder width.

**Indications**

Non-surgical clavicle fractures.

**Contraindications**

Do not use the product if the diagnosis has not been confirmed.

Do not use in the event of known allergy to any of the components.

Do not apply the product in direct contact with broken skin.

Do not use in the event of complex or surgical shoulder fractures or rib cage fractures.

Do not use in case of severe vascular conditions of the upper limb (thrombosis, lymphedema, arterial compression...).

**Precautions**

Verify the product's integrity before every use.

Do not use the product if it is damaged.

Choose the appropriate size to fit the patient, referring to the size chart.

Wear the device above thin clothing.

It is recommended that a healthcare professional supervises the first application.

Strictly comply with your healthcare professional's prescription and recommendations for use.

For hygiene, security and performance reasons, do not re-use the product for another patient.

In the event of discomfort, significant hindrance, pain, variation in limb volume, abnormal sensations or change in colour of the extremities, remove the device and consult a healthcare professional.

It is recommended to adequately tighten the device to achieve support/immobilisation without restricting blood circulation.

Do not use the device in case of application of certain products on the skin (creams, ointments, oils, gels, patches...).

Do not wear the product when driving a vehicle.

The use of a medical device by a child, should be done under the supervision of an adult or a healthcare professional.

This product is intended for the treatment of a given condition. Its duration of use is limited to this treatment only.

**Undesirable side-effect**

This device can cause skin reactions (redness, itching, burns, blisters, etc.) or wounds of various degrees of severity.

Possible risks of clavicle fracture complications.

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

**Instructions for use/Application**

The help of another person is required to fit this device.

1. Remove the self-fastening Y-tabs from the edges of the straps, trim the straps according to the patient's shoulder width and along the corresponding visual marker. Then, reposition the self-fastening Y-tabs on the edges of the straps.
2. Centre the dorsal anchoring pad. Adjust the height if needed. Position the straps around the collar bones and under the armpits.
3. Pass the straps through the buckles, adjust the tightening and fasten the straps via the self-fastening tabs.

**Care/Maintenance**

Product can be washed in accordance with the instructions shown on this leaflet and on the label. Hand

washing. Remove the self-fastening tabs before washing. Do not dry clean. Do not tumble-dry. Do not iron. Do not use detergents, fabric softeners or aggressive products (products containing chlorine). Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.).

**Storage**

Store at room temperature, preferably in the original packaging.

**Disposal**

Dispose of in accordance with local regulations.

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© Thuasne - 6172201 (2022-12)
Thuasne UK Ltd