
THUASNE®


Dynabelt

fr	 Ceinture de soutien abdominal.....	5
en	 Abdominal support belt.....	6
de	 Bauch-Stabilisierungsorthese	8
nl	 Abdominale brace.....	10
it	 Cintura di sostegno addominale	12
es	 Faja de sujeción abdominal.....	14
pt	 Cinta de suporte abdominal	16
da	 Mavestøttebælte	18
fi	 Vatsan tukivyö.....	20
sv	 Stödjande maggördel.....	22
el	 Ζώνη υποστήριξης της κοιλιακής χώρας.....	23
cs	 Břišní podpurný pás.....	25
pl	 Pas podtrzymujący mięśnie brzucha.....	27
lv	 Vēdera daļas atbalsta josta	29
lt	 Palaikomasis pilvo diržas.....	31
et	 Kõhupiirkonna tugikorsett.....	33
sl	 Oporni pas za trebušni predel	34
sk	 Podporný brušný pás.....	36
hu	 Haskötő.....	38
bg	 Поддържащ колан за коремната област	40
ro	 Centură de susținere abdominală	42
ru	 Опорный брюшной корсет	44
hr	 Abdominalni potporni pojasi	46
zh	 腹部稳定护带	48
ar	 حزام لدعم البطن.....	51


Instructions d'utilisation - Instructions for use - Gebrauchsanweisung - Gebruiksaanwijzingen - Istruzioni per l'uso - Instrucciones de uso - Instruções de utilização - Anvisninger vedrørende brug - Käyttöohjeet - Råd vid användning - Οδηγίες χρήσης - Pokyny k používání - Szczegóły użytkowania - Lietošanas instrukcija - Naudojimo instrukcijos - Kasutusjuhised - Navodila za uporabo - Pokyny na používanie - Használati utasítás - Инструкции за употреба - Instrucțiuni de utilizare - Инструкция по эксплуатации - Upute za uporabu - 使用指南 - إرشادات الاستعمال

(A)


(B)


(C)


(D)

(E)

	1	60 - 80 cm
	2	80 - 100 cm
	3	100 - 130 cm*
	4	130 - 160 cm*

fr	Tour de taille - Hauteurs : 18 cm/25 cm/33 cm *Pour la France uniquement, sur-mesure patient.
en	Waist circumference - Heights: 18 cm/25 cm/33 cm *For France only, made to measure patient.
de	Taillenumfang - Höhen: 18 cm/25 cm/33 cm *Nur für Frankreich, maßgeschneidert.
nl	Tailleomvang - Hoogtes: 18 cm / 25 cm / 33 cm *Enkel voor Frankrijk, op maat van de patiënt.
it	Circonferenza vita - Altezze: 18 cm/25 cm/33 cm *Solo per la Francia, su misura per il paziente.
es	Contorno de cintura - Alturas: 18 cm/25 cm/33 cm *Solo para Francia, a medida del paciente.
pt	Perímetro de cintura - Alturas: 18 cm/25 cm/33 cm *Apenas para a França, à medida do paciente.
da	Taljemål - Højde: 18 cm/25 cm/33 cm *Kun for Frankrig, efter patientens mål.
fi	Vyötärönympärys - Korkeudet: 18 cm/25 cm/33 cm *Ainoastaan Ranskassa, potilaan mittojen mukaan.
sv	Midjemått - Höjd: 18 cm/25 cm/33 cm *Endast för Frankrike, patientspecifik.
el	Περίμετρος μέσης - Ύψη: 18 cm/25 cm/33 cm *Μόνο για Γαλλία, κατά παραγγελία ασθενούς.
cs	Obvod pasu - Výšky: 18 cm / 25 cm / 33 cm *Výhradně pro Francii, na míru pacienta.
pl	Obwód w talii - Wysokości: 18 cm/ 25 cm/ 33 cm *Tylko we Francji produkt jest wykonywany na miarę dla pacjenta.
lv	Vidukļa apkārtmērs - Augstumi: 18 cm/25 cm/33 cm *Tikai attiecībā uz Franciju - pacientam pielāgots izstrādājums.
lt	Juosmens apimtis - Aukštis: 18 cm / 25 cm / 33 cm *Tik Prancūzijai, pagal paciento užsakymą.
et	Vööümberrmõõt - Kõrgused: 18 cm/25 cm/33 cm *Ainult Prantsusmaal, patsiendile kohandatud.
sl	Obseg pasu - Višina: 18 cm/25 cm/33 cm *Samo za Francijo, po meri pacienta.
sk	Obvod pása - Výšky: 18 cm/25 cm/33 cm *Len pre Francúzsko, vyrobené na mieru pacienta.
hu	Derékbőség - Magasság: 18 cm/25 cm/33 cm *Csak Franciaországban, a páciens egyéni méretére szabva.
bg	Обиколка на талията - Височини: 18 cm/25 cm/33 cm *Само за Франция, изработена по мярка за пациента.
ro	Circonfenrința taliei - Înălțimi: 18 cm/25 cm/33 cm *Doar pentru Franța, făcut pe măsura pacientului.
ru	Обхват талии - Высота: 18 см/25 см/33 см *Только для Франции, изготавливается по меркам пациента.
hr	Opseg struka - Visine: 18 cm/25 cm/33 cm *Samo u Francuskoj, po mjeri bolesnika.
zh	腰围 - 高度 : 18厘米/25厘米/33厘米 *仅限于法国, 为患者定制。
ar	محيط الخصر - الارتفاعات: 18 سم / 25 سم / 33 سم *موجه لفرنسا فقط، حسب مقياس المريض

fr	Compression
en	Compression
de	Kompression
nl	Compressie
it	Compressione
es	Compresión
pt	Compressão
da	Kompression
fi	Kompressio
sv	Kompression
el	Συμπίεση
cs	Komprese
pl	Kompresja
lv	Kompresija
lt	Kompresija
et	Kompressioon
sl	Kompresija
sk	Kompresia
hu	Szorítás
bg	Компресия
ro	Compresie
ru	Компрессия
hr	Kompresija
zh	压紧
ar	الضغط

Maintien
Support
Halt
Fixatie
Sostegno
Sujeción
Manutenção
Støtte
Tuki
Stöd
Συγκράτηση
Stabilizace
Podtrzymywanie
Atbalsts
Palaikymas
Hooldus
Oprijem
Podpora
Támasz
Поддръжка
Susținere
Стабилизация
Potpورا
稳定
الدعم

CEINTURE DE SOUTIEN ABDOMINAL

Description/Destination

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ABDOMINAL SUPPORT BELT

Description/Destination

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table*.

Composition

Fastening "loop knit" elastic fabric. ❶

Self-fastening with hand loop. ❷

Rigid, conformable and removable dorsal stays. ③

Flexible and removable abdominal stays. ④

Hypogastric strap. ⑤

Textile components: polyurethane - elastane - polyamide - polyester - cotton.

Non-textile components: brushed tempered steel - polyoxymethylene - polyamide.

Properties/Mode of action

The belt ensures strong support and compression suitable for the abdomen thanks to the elastic fabric and the flexible stays. The removable hypogastric strap provides extra support around the abdomen if needed.

Indications

Abdominal muscle weakness.

Postoperative conditions: laparotomy laparoscopy.

Postpartum (after a cesarean).

Orthostatic hypotension.

Contraindications

Do not apply the product in direct contact with broken skin.

Do not use in the event of known allergy to any of the components.

Do not use on pregnant women.

Do not use in the event of a hiatus hernia.

Do not use in the event of bone cancer with spinal metastases.

Do not use in the event of circulatory, pulmonary or cardiovascular problems in patients for whom an increase in arterial blood pressure is not recommended.

Surgeons usually do not recommend any kind of compression for a period of time after cardiac or thoracic surgery for people at risks.

Do not use in patients fed by a gastrostomy tube or with aortic aneurysm.

Precautions

Verify the product's integrity before every use.

Do not use the device if it is damaged.

Choose the appropriate size to fit the patient, referring to the size chart.

Strictly comply with your healthcare professional's prescription and recommendations for use.

It is recommended that a healthcare professional adjust the stays to the shape of the patient's back before the belt is worn for the first time.

It is recommended that a healthcare professional supervises the first application.

Consult a healthcare professional before use in case of pelvic floor deficiency.

Do not wear the product in a medical imaging device.

In the event of discomfort, significant hindrance, pain, abnormal sensations, remove the device and consult a healthcare professional.

This product is intended for the treatment of a given condition. Its duration of use is limited to this treatment only.

For hygiene and performance reasons, do not re-use the product for another patient.

It is recommended to adequately tighten the device to achieve support without excessive compression.

Do not use the device in case of application of certain products on the skin (creams, ointments, oils, gels, patches...).

The use of a medical device by a child, should be done under the supervision of an adult or a healthcare professional.

Undesirable side-effect

This device can cause skin reactions (redness, itching, burns, blisters, etc.) or wounds of various degrees of severity.

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Application

Position the belt around the waist, using the hand-loops. Ⓐ

The belt must be centred on the spine: the dorsal stays should be positioned on either side of the spine.

The bottom of the belt must stop at the top of the crease between the buttocks.

Hold the end without a hook on the abdomen. Ⓑ

Hold the other end and place the hook onto the fastening "loop knit" elastic fabric. Adjust to the desired level of tightness. Ⓒ

The level of compression can be adjusted as required throughout the day.

If necessary, position the hypogastric strap on the ventral part of the Dynabelt Ⓓ, then pull on its two ends and secure them once the desired level of support is obtained Ⓔ.

Care/Maintenance

Close the self-fastening tabs before washing. Machine washable at 30°C. Do not machine wash the device more than 100 times. Do not dry clean. Do not use detergents, fabric softeners or aggressive products (products containing chlorine). Do not tumble-dry. Do not iron. Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.). Dry flat.

Storage

Store at room temperature, preferably in the original packaging.

Disposal

Dispose of in accordance with local regulations.

* Not applicable to made-to-measure medical devices, where measurements will be taken individually for each patient.

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Mantenga el segundo extremo y coloque el corchete en el tejido elástico "rizo" adherente. Ajuste el nivel de apriete deseado. © En cualquier momento del día, podrá ajustar el nivel de compresión en función de sus necesidades.
THUASNE UK Ltd
Thuasne SAS - SIREN/RCS Nanterre 542 091 186