



THUASNE®

MEDICAL COMPRESSION

INTERNATIONAL

2026

VENOFLEX

Biflex & Kit Biflex



Editorial

Dear Partners,

Thuasne®'s mission is to empower each patient and caregiver in their efforts to restore mobility and improve quality of life through non-invasive care. By placing prevention and independence at the heart of our approach, we offer a holistic vision of health, where each individual becomes an active player in the healing process.

A family-owned business founded in 1847, we are now an integrated French mid-sized industrial company, controlling the entire value chain – from thread to delivery. With a strong territorial network, we proudly uphold an industrial model that creates jobs, drives innovation, and develops local expertise.

Because we have a demanding vision of healthcare – high-quality healthcare, without compromise – we are committed to working alongside you to support an intelligent, sustainable, and responsible system.

Looking ahead to 2026, our ambition remains unchanged: to provide your patients with innovative solutions with proven therapeutic effectiveness, and to demonstrate their value through clinical, biomechanical, and health economics studies. To this day, more than 20 studies are underway to assess the medical benefits of our solutions.

We are also actively advancing our digital transformation. We have launched a measurement application designed to deliver devices tailored to each patient. At the same time, integrating artificial intelligence into our tools allows us to focus our efforts on what truly matters to you – where our expertise provides real added value.

Being a healthcare player also means being responsible. Through our initiative “Move for Tomorrow!”, we are strengthening our CSR commitments: removing or replacing plastics in our packagings, reducing energy consumption at our industrial sites by 5.4%*, supporting the Mont-Blanc accessibility tour... These are just a few concrete steps toward building more sustainable healthcare and a better future.

In 2025, you discovered Biflex Self Adjust and 2 new colors (Matcha & Pivoine) for Venoflex Fast Lin, the result of the expertise of our R&D, marketing, design, and production teams. In 2026, we will continue to expand our range of products, services, and training to better support your needs. The year will be marked by the launch of new solutions that we are eager to share with you.

We warmly thank you for your continued trust and wish you a wonderful year in 2026.

Delphine Hanton
Chief Executive Officer



Elizabeth Ducottet
Chair Woman



175 years of the Thuasne® Group

The know-how of the Thuasne® Group is rooted in Saint-Étienne, having imagined, developed and produced medical devices for over 175 years, allowing everyone to take charge of their own health.

As a result of the shared history of 6 generations of entrepreneurs, Thuasne® offers specific, simple, adapted and innovative health solutions supporting mobility in over 100 countries. We offer our expertise to patients in the fields of orthopaedics (belts, orthoses), medical compression (phlebology, lymphology, healing), femininity (external breast prostheses, underwear) and home care.

"Made in Thuasne®": the ability to control the whole value chain

Because Thuasne® pursues the objective of offering patients and healthcare professionals effective, quality products, it can guarantee that its products are designed, produced and distributed using technologies controlled by the Group, after extensive accumulated experience.

Today, Thuasne® is...



2,500
employees



Revenue of
€290M
in 2024



50%
of revenue earned
internationally



Commercial presence in
110 countries



17
industrial sites in Europe,
North America, and Africa



23,000
B2B customers



10,000
orders handled
per day worldwide



5
international patents
registered each year



21 millions
km of knitted textiles



5,300
km of material produced per year



A Thuasne® brace is sold every
6 seconds
in the world



4
R&D centers

Contents

Ceap classification* and Thuasne® solutions.....	/ p. 5
Taking measurements.....	/ p. 6
Range recap VENOFLEX WOMEN	/ p. 10
Range recap VENOFLEX MEN	/ p. 14
Thuasne® recommends.....	/ p. 16
Measuring guide.....	/ p. 17
How do you deliver the right solution?.....	/ p. 18

VENOFLEX

Collection WOMEN.....	/ p. 24	Maternity	/ p. 76
Collection MEN.....	/ p. 32	Specific UNISEX.....	/ p. 84
Discreet WOMEN.....	/ p. 36	Original WOMEN.....	/ p. 90
Discreet MEN.....	/ p. 48	Original MATERNITY.....	/ p. 96
Sensitive WOMEN.....	/ p. 52	Original MEN.....	/ p. 100
Sensitive MEN.....	/ p. 62		

Accessories & adapted products VENOFLEX	/ p. 104
Bandages.....	/ p. 106
Other Thuasne® expertise.....	/ p. 115

Availability of these products might vary from a given country or region to another.
All the medical devices class I, II, mentioned on this document are CE marked according to the European Regulation 2017/745 on medical devices, unless specifically identified as "NOT CE MARKED".
Please read carefully the instructions for use of the products.

 Thuasne - 120, rue Marius Auphan 92300 - Levallois Perret - France.
www.thuasne.com

Last revision date: 02/2026



CEAP CLASSIFICATION* AND THUASNE® SOLUTIONS

CLASSIFICATION OF CHRONIC VENOUS DISEASES

	C0	C1	C2	C3	C4a	C4b	C4c	C5	C6
Indications	No visible signs Pain, heaviness, cramps	Telangiectasia, reticular veins < 3 mm	Varicose veins > 3 mm	Oedema	Pigmentation, venous eczema	Lipodermato- sclerosis	Corona phlebectatica	Healed ulcer	Open ulcer
RECOMMENDED ASQUAL NORM	1 10 – 15 mmHg	1 10 – 15 mmHg 2 15 – 20 mmHg	2 15 – 20 mmHg 3 20 – 36 mmHg	3 20 – 36 mmHg	3 20 – 36 mmHg	3 20 – 36 mmHg	3 20 – 36 mmHg	3 20 – 36 mmHg 4 > 36 mmHg	4 > 36 mmHg
RECOMMENDED RAL NORM		1 18 – 21 mmHg	1 18 – 21 mmHg 2 23 – 32 mmHg	2 23 – 32 mmHg	2 23 – 32 mmHg	2 23 – 32 mmHg	2 23 – 32 mmHg	2 23 – 32 mmHg 3 34 – 46 mmHg 4 > 49 mmHg	3 34 – 46 mmHg 4 > 49 mmHg

THUASNE® SOLUTIONS


ASQUAL	1 10 – 15 mmHg	2 15 – 20 mmHg		3 20 – 36 mmHg
	VENOFLEX <i>Secret</i>	VENOFLEX <i>Fantaisie</i> • OGEE	VENOFLEX <i>Fast</i> LIN	VENOFLEX <i>Fast</i> LAINE
	VENOFLEX <i>City Confort</i> COTON	VENOFLEX <i>Incognito</i> ABSOLU	VENOFLEX <i>Fast</i> LAINE	VENOFLEX <i>Fast</i> AIR
	VENOFLEX <i>Elégance</i>	VENOFLEX <i>Incognito</i> ABSOLU • SWISS DOTS	VENOFLEX <i>Fast</i> AIR	VENOFLEX <i>Fast</i> COTON
	VENOFLEX <i>Clinic</i>	VENOFLEX <i>Kokoon</i> ABSOLU	VENOFLEX <i>City Confort</i> COTON	VENOFLEX <i>Kokoon</i> ABSOLU
		VENOFLEX <i>Kokoon</i> ABSOLU • LOZENGES	VENOFLEX <i>Secret</i>	VENOFLEX <i>City Confort</i> COTON
		VENOFLEX <i>Kokoon</i> ABSOLU • MATERNITY	VENOFLEX <i>Secret</i> • MATERNITY	VENOFLEX <i>Secret</i>
		VENOFLEX <i>City Confort</i> FIL D'ÉCOSSE	VENOFLEX <i>Elégance</i>	VENOFLEX <i>Elégance</i>
		VENOFLEX <i>Simply</i> COTON FIN	VENOFLEX <i>Soft & Care</i>	VENOFLEX <i>Soft & Care</i>
		VENOFLEX <i>Fast</i> COTON	VENOFLEX <i>Clinic</i>	

RAL	1 18 – 21 mmHg	2 23 – 32 mmHg	
	VENOFLEX <i>Micro</i> RAL	VENOFLEX <i>Micro</i> RAL	from C3
		VENOFLEX <i>Men</i> RAL	from C4a
		VENOFLEX <i>Elegance</i> RAL	from C5
			from C6

BANDAGE BIFLEX 16+ PRATIC BANDAGE BIFLEXIDEAL
COTTON BANDAGE**
VENOFLEX MADE-TO-MEASURE PATIENT layering of 15 – 20 mmHg VENOFLEX stockings
Kit Biflex

* Source : Eklöf B : Revision of the CEAP classification for chronic venous disorders: consensus statement ; J Vasc Surg 2004.

** To be used in addition to a compressive layer.



- cB Ankle circumference
(at narrowest point)
- cC Calf circumference (at widest point)
- cG Thigh circumference (at widest point)
- cH Hip circumference
- cT Waist circumference
- lD Floor-to-knee length
(measurement taken below
the popliteal fossa, posterior aspect)
- lK Floor-to-crotch length
(internal face of the leg)



www.asqual.com

ASQUAL

- CLASS ① 10 – 15 mmHg
- CLASS ② 15 – 20 mmHg
- CLASS ③ 20 – 36 mmHg

SOCKS (AD)



measurements in cm			cB	cC-	cC	cC+
WOMEN	NORMAL ℓD < 40 LONG ℓD > 40	0	16 – 19	21 – 28	28 – 33	33 – 42
		1	19 – 21	23 – 30	30 – 35	35 – 46
		2	21 – 24	26 – 32	32 – 38	38 – 48
		3	24 – 26	29 – 34	34 – 42	42 – 52
		4	26 – 29	32 – 37	37 – 46	46 – 54
		5	29 – 31	36 – 43	43 – 49	49 – 56
MEN	NORMAL ℓD < 44 LONG ℓD > 44	1	19 – 21	23 – 30	30 – 35	35 – 46
		2	21 – 24	26 – 32	32 – 38	38 – 48
		3	24 – 26	29 – 34	34 – 42	42 – 52
		4	26 – 29	32 – 37	37 – 46	46 – 54
		5	29 – 31	36 – 43	43 – 49	49 – 56

STOCKINGS (AG)



		measurements in cm	cB	cG-	cG	cG+
WOMEN	NORMAL ℓ _K < 74 LONG ℓ _K > 74	0	16 – 19	35 – 42	42 – 52	52 – 57
		1	19 – 21	38 – 45	45 – 54	54 – 59
		2	21 – 24	43 – 50	50 – 59	59 – 66
		3	24 – 26	48 – 56	56 – 66	66 – 70
		4	26 – 29	51 – 58	58 – 69	69 – 74
		5	29 – 31	55 – 63	63 – 74	74 – 78
MEN	NORMAL ℓ _K < 78 LONG ℓ _K > 78	1	19 – 21	-	45 – 54	-
		2	21 – 24	-	50 – 59	-
		3	24 – 26	-	56 – 66	-
		4	26 – 29	-	58 – 69	-

WOMEN

To fit shoe size: 36 – 42 (UK 3 – 8).
Size 0 available in Venoflex Kokoon Absolu only.
Socks +/- , stockings +/- and size 5 available
in Venoflex Kokoon Absolu and Venoflex Secret only.

MEN

To fit shoe size: 38 – 46 (UK 5 – 11).
Long feet socks (EU 45 – 48) available in Venoflex City
Confort Fil d'Ecosse and Venoflex Elegance.
Socks +/- available in Venoflex City Confort Fil d'Ecosse
and Venoflex Elegance.
Size 5 available in Venoflex City Confort Fil d'Ecosse.



		measurements in cm		cB	cG	cH max.	cT max.
TIGHTS (AT)	WOMEN	NORMAL $\ell_K < 74$ LONG $\ell_K > 74$	0	16 - 19	42 - 52	120	90
			1	19 - 21	45 - 54	128	93
			2	21 - 24	50 - 59	133	96
			3	24 - 26	56 - 66	145	104
			4	26 - 29	58 - 69	154	110
			5	29 - 31	63 - 74	162	113
	MEN	NORMAL $\ell_K < 78$ LONG $\ell_K > 78$	1	19 - 21	45 - 54	120	100
			2	21 - 24	50 - 59	135	115
			3	24 - 26	56 - 66	135	115
			4	26 - 29	58 - 69	155	135



		measurements in cm		cB	cG	cH max.	cT max.
MATERNITY (AT)	WOMEN	NORMAL $\ell_K < 74$ LONG $\ell_K > 74$	0	16 - 19	42 - 52	136	105
			1	19 - 21	45 - 54	136	105
			2	21 - 24	50 - 59	145	115
			3	24 - 26	56 - 66	158	125
			4	26 - 29	58 - 69	167	130
			5	29 - 31	63 - 74	176	135

WOMEN
To fit shoe size: 36 - 42 (UK 3 - 8).
Size 0 available in Venoflex Kokoon Absolu only.

MEN
To fit shoe size: 38 - 46 (UK 5 - 11).



		measurements in cm		cB	cC
ANTI-STASIS SOCKS	UNISEX	NORMAL $\ell_D < 40$ LONG $\ell_D > 40$	1	17 - 19	30 - 35
			2	19 - 22	32 - 38
			3	22 - 25	34 - 42
			4	25 - 29	37 - 46
			5	29 - 31	43 - 49



		measurements in cm		cB	cG
ANTI-STASIS STOCKINGS	UNISEX	NORMAL $\ell_K < 74$ LONG $\ell_K > 74$	1	17 - 19	45 - 54
			2	19 - 22	50 - 59
			3	22 - 25	56 - 66
			4	25 - 29	58 - 69
			5*	29 - 31	63 - 74

*15 - 20 mmHg only.
To fit shoe size 36 - 46 (UK 3 - 11).


VENOFLEX
Fantaisie · OGEE



		measurements in cm		cB	cG
STOCKINGS (AG)	WOMEN	NORMAL $\ell_K < 74$ LONG $\ell_K > 74$	1	19 - 21	50 - 59
			2	21 - 24	55 - 64
			3	24 - 26	61 - 70



		measurements in cm		cB	cG	cH max.	cT max.
TIGHTS (AT)	WOMEN	NORMAL $\ell_K < 74$ LONG $\ell_K > 74$	1	19 - 21	50 - 59	128	93
			2	21 - 24	55 - 64	133	96
			3	24 - 26	61 - 70	145	104



- cA Circumference of the head of the metatarsal bones
- cY Instep circumference
- cB Ankle circumference (at narrowest point)
- cB1 Circumference at the start of the calf muscles
- cC Calf circumference (at widest point)
- cD Circumference below the knee
- cE Knee circumference
- cF Circumference at the equidistant point between G and E
- cG Thigh circumference (at widest point)
- cH Hip circumference
- cT Waist circumference
- lD Floor-to-knee length (measurement taken below the popliteal fossa, posterior aspect)
- lK Floor-to-crotch length (internal face of the leg)
- lG Floor-to-gluteal fold length (external face of the leg)



- RAL** CLASS ① 18 – 21 mmHg
- CLASS ② 23 – 32 mmHg

SOCKS (AD)



measurements in cm			cA	cY	cB	cB1	cC	cD
UNISEX	SHORT D 34-39 LONG D 39-45	1	17 - 22	26 - 31	18 - 20	23 - 27	28 - 34	27 - 33
		2	19 - 24	28 - 33	20 - 22	24 - 29	30 - 37	29 - 36
		3	21 - 26	29 - 35	22 - 24	26 - 32	33 - 40	32 - 39
		4	23 - 29	31 - 37	24 - 26	29 - 35	35 - 43	34 - 42
		5	25 - 32	32 - 38	26 - 28	31 - 37	37 - 46	36 - 45
		6	27 - 34	33 - 40	28 - 30	33 - 39	39 - 49	38 - 48

STOCKINGS (AG)



measurements in cm			cA	cY	cB	cB1	cC	cD	cE	cF	cG	cG+
UNISEX	SHORT LG 60-69 NORMAL LG 70-80 LONG LG 81-90	1	17 - 22	26 - 31	18 – 20	23 - 27	28 - 34	27 - 33	30 - 37	39 - 52	43 - 48	49 - 57
		2	19 - 24	28 - 33	20 - 22	24 - 29	30 - 37	29 - 36	33 - 40	41 - 56	45 - 52	53 - 62
		3	21 - 26	29 - 35	22 - 24	26 - 32	33 - 40	32 - 39	35 - 43	44 - 60	49 - 56	57 - 67
		4	23 - 29	31 - 37	24 - 26	29 - 35	35 - 43	34 - 42	37 - 45	47 - 65	53 - 60	61 - 72
		5	25 - 32	32 - 38	26 - 28	31 - 37	37 - 46	36 - 45	39 - 48	50 - 69	56 - 64	65 - 77
		6	27 - 34	33 - 40	28 - 30	33 - 39	39 - 49	38 - 48	41 - 51	53 - 73	60 - 68	69 - 81

TIGHTS (AT)



measurements in cm			cA	cY	cB	cB1	cC	cD	cE	cF	cG	cH max.	cT max.
UNISEX	NORMAL ℓ_G 65-75 LONG ℓ_G 75-85	1	17 - 22	26 - 31	18 - 20	23 - 27	28 - 34	27 - 33	30 - 37	39 - 52	43 - 57	85 - 105	59 - 70
		2	19 - 24	28 - 33	20 - 22	24 - 29	30 - 37	29 - 36	33 - 40	41 - 56	45 - 62	90 - 110	66 - 80
		3	21 - 26	29 - 35	22 - 24	26 - 32	33 - 40	32 - 39	35 - 43	44 - 60	49 - 67	95 - 115	76 - 85
		4	23 - 29	31 - 37	24 - 26	29 - 35	35 - 43	34 - 42	37 - 45	47 - 65	53 - 72	100 - 120	82 - 90
		5	25 - 32	32 - 38	26 - 28	31 - 37	37 - 46	36 - 45	39 - 48	50 - 69	56 - 77	105 - 125	88 - 95
		6	27 - 34	33 - 40	28 - 30	33 - 39	39 - 49	38 - 48	41 - 51	53 - 73	60 - 81	110 - 130	92 - 105

MATERNITY (AT)



measurements in cm			cA	cY	cB	cB1	cC	cD	cE	cF	cG	cH max.	cT max.
WOMEN	NORMAL ℓG 65-75	1	17 - 22	26 - 31	18 - 20	23 - 27	28 - 34	27 - 33	30 - 37	39 - 52	43 - 57	85 - 105	59 - 70
		2	19 - 24	28 - 33	20 - 22	24 - 29	30 - 37	29 - 36	33 - 40	41 - 56	45 - 62	90 - 110	66 - 80
		3	21 - 26	29 - 35	22 - 24	26 - 32	33 - 40	32 - 39	35 - 43	44 - 60	49 - 67	95 - 115	76 - 85
	LONG ℓG 75-85	4	23 - 29	31 - 37	24 - 26	29 - 35	35 - 43	34 - 42	37 - 45	47 - 65	53 - 72	100 - 120	82 - 90
		5	25 - 32	32 - 38	26 - 28	31 - 37	37 - 46	36 - 45	39 - 48	50 - 69	56 - 77	105 - 125	88 - 95
		6	27 - 34	33 - 40	28 - 30	33 - 39	39 - 49	38 - 48	41 - 51	53 - 73	60 - 81	110 - 130	92 - 105

VENOFLEX
Men

SOCKS (AD)



measurements in cm			cA	cY	cB	cB1	cC	cD
MEN	SHORT ℓ _D < 44,5 LONG ℓ _D ≥ 44,5	0	16 - 21	24 - 30	17 - 19	22 - 25	27 - 35	26 - 34
		1	18 - 23	27 - 32	19 - 21	23 - 28	30 - 35	29 - 34
		2	20 - 26	28 - 35	21 - 24	25 - 32	32 - 38	31 - 37
		3	23 - 29	31 - 37	24 - 26	29 - 35	34 - 42	33 - 41
		4	25 - 33	32 - 39	26 - 29	31 - 38	37 - 46	36 - 45
		5	28 - 35	34 - 41	29 - 31	34 - 40	43 - 49	42 - 48

WOMEN

ASQUAL

Collection

Discreet

Sensitive

Maternity

Specific

Original

			CLASS	SOCKS	SOCKS CALF
VENOFLEX Fantaisie · OGEE	p. 26	②			
VENOFLEX Incognito ABSOLU · SWISS DOTS	p. 28	②	black swiss dots		
VENOFLEX Kokoon ABSOLU · LOZENGES	p. 30	②	black diamonds		
VENOFLEX Incognito ABSOLU	p. 38	②	natural nude golden bronzed dark brown black		
VENOFLEX Kokoon ABSOLU	p. 42	②	navy honey natural beige black natural beige black		
VENOFLEX Fast LIN	p. 54	②	azure blue black strips natural strips		
VENOFLEX Fast COTON	p. 56	②	grey irish marl black irish		
VENOFLEX Fast LAINE	p. 58	②	natural twists anthracite twists black twists		
VENOFLEX Simply COTON FIN	p. 60	②	natural beige grey violet powder pink cream black		
VENOFLEX Kokoon ABSOLU · MATERNITY	p. 78	②			
VENOFLEX Soft & Care	p. 86	②	black		
VENOFLEX Clinic	p. 88	①			
VENOFLEX Secret	p. 92	②	golden black navy natural beige dune golden golden dark brown NEW black tan beige dark brown NEW black		
VENOFLEX Secret OPAQUE	p. 94	②	golden black		
VENOFLEX Secret · MATERNITY	p. 98	②	navy golden black		

RAL



GÜTEZEICHEN



MEDIZINISCHE
KOMPRESSIONS-
STRÜMPFE

WOMEN

RAL

CLASS



SOCKS



STOCKINGS



STOCKINGS
⊖ ⊕



TIGHTS

VENOFLEX
Micro

p. 44

1

light caramel tan beige black

light caramel tan beige black

light caramel tan beige black

light caramel tan beige black

2

light caramel tan beige black

light caramel tan beige black

light caramel tan beige black

light caramel tan beige black

VENOFLEX
Elegance

p. 46

2

light caramel tan beige diamond
marine black

light caramel tan beige diamond
marine black

light caramel tan beige diamond
marine black

light caramel tan beige diamond
marine black

Discreet

CLASS



MATERNITY TIGHTS
& FULL HIPS

VENOFLEX
Micro

p. 80

1

light caramel tan beige black

2

light caramel tan beige black

VENOFLEX
Elegance

p. 82

2

light caramel tan beige diamond
marine black

Maternity

Define your patients' needs for better compliance

COMPRESSION
CLASS

ASQUAL

Class
1
10 - 15 mmHg
2
15 - 20 mmHg
3
20 - 36 mmHg

RAL

Class
1
18 - 21 mmHg
2
23 - 32 mmHg

FORMAT



SOCKS



STOCKINGS



TIGHTS



MATERNITY
TIGHTS



Collection

Sensitive

Specific

Original

MEN

ASQUAL

			CLASS	SOCKS	SOCKS CALF
	VENOFLEX <i>City Confort</i> FIL D'ÉCOSSE	p. 34	②	beige dakar black navy blue	black
	VENOFLEX <i>Fast</i> AIR	p. 64	② ③	navy anthracite black black	
	VENOFLEX <i>Fast</i> LIN	p. 66	②	azure blue greige natural black	
	VENOFLEX <i>Fast</i> COTON	p. 68	② ③	slate ribbed black ribbed black ribbed	
	VENOFLEX <i>Fast</i> LAINE	p. 70	② ③	anthracite black black	
	VENOFLEX <i>City Confort</i> COTON	p. 72	① ② ③	black granite black navy blue black	
	VENOFLEX <i>Soft & Care</i>	p. 86	② ③	black black	
	VENOFLEX <i>Clinic</i>	p. 88	① ②		
	VENOFLEX <i>Elégance</i>	p. 102	① ② ③	black beige sand navy grey black beige sand black	black

Discreet

Sensitive

MEN

RAL

			CLASS	SOCKS	STOCKINGS	STOCKINGS CALF
	VENOFLEX <i>Elegance</i>	p. 50	②	light caramel tan beige diamond marine black	light caramel tan beige diamond marine black	light caramel tan beige diamond marine black
	VENOFLEX <i>Men</i>	p. 74	②	black navy blue		

 OPEN-TOE SOCKS	 OPEN-TOE STOCKINGS	 TIGHTS
 black	 black	
		 black
 black	 black	
	 white	
 white	 white	
 beige sand	 black	
 beige sand	 black	

 TIGHTS
     

From caring for your patients to delivering the VENOFLEX range

A GOOD
MEASUREMENT

=
*effective
compression*

Measurements should be taken in the morning,
to avoid swollen legs, which could distort the
measurements.

A PRODUCT
ADAPTED TO THE PATIENT

=
*better
compliance*

Once you have reviewed the prescription, here
are the steps you need to take to best meet your
patients' needs:



5

**KEY
STEPS:**

1 *Patient consultation*

background, symptoms, frequency, duration, severity, prescription analysis.

2 *Measurements*

of both legs, even at each renewal.

3 *Questions*

about their lifestyle, their activities, their tastes, their sensitivity to materials, their sensitivity to heat...

4 *Product and size selection*

Beware of **seasonal changes**

Remember to offer a wide range of venous compression products: opacity or transparency, light fibers (linen, cotton) or more wintry fibers (wool), open or closed feet, various colors... so that your patients can wear compression in all seasons, to avoid the risk of venous insufficiency symptoms reappearing, and promote compliance.

5 *Donning, fitting and care advice*

Beware of **wear and tear**

Washing, rubbing, aging fibers and their elasticity mean that compression weakens. Make your patients aware of the need to pay attention to their compression stockings to maintain the product's effectiveness.

Beware of **darning**


Raise your patients' awareness: if socks have holes in them, darning can **compromise the product's efficacy**.

PROMOTE COMPLIANCE

- Washable at  40
- Do not tumble dry
- Do not add fabric softener
- Remove jewelry before donning




A donning sock is provided in the boxes
of Class **3** products.



How to measure your patients properly?

IDEAL CONDITIONS FOR MEASUREMENTS⁽¹⁾:

- In a confidential area.
- Preferably in the morning when the legs are not swollen.
- On both bare legs, feet flat on the floor.
- In case of mild to moderate oedema: elevate the legs for 10 to 20 min.
- Define the measurements below, then refer to the sizing chart on the product page.




Products under the ASQUAL norm only require the measurements with this symbol: ➤.

- cA Circumference of the head of the metatarsal bones
- cY Instep circumference
- cB Ankle circumference (at narrowest point)
- cB1 Circumference at the start of the calf muscles
- cC Calf circumference (at widest point)
- cD Circumference below the knee
- cE Knee circumference
- cF Circumference at the equidistant point between G and E
- cG Thigh circumference (at widest point)
- cH Hip circumference
- cT Waist circumference
- lD Floor-to-knee length (measurement taken below the popliteal fossa, posterior aspect)
- lK Floor-to-crotch length (internal face of the leg)
- lG Floor-to-gluteal fold length (external face of the leg)



Renewal⁽¹⁾:

Repeat measurements on both bare legs.



(1) Fiche compression A4 produit, pour l'information du patient en traitement, proposé par la SFMV, juillet 2022. https://www.portailvasculaire.fr/sites/default/files/imce/Fiches%20pratiques/Fiches_patients_compression/0-FICHE_1COMPRESSION_PATIENT.pdf

Use and care of VENOFLEX socks, stockings and tights

WARNING ABOUT MORPHOLOGICAL CHANGES

All morphological changes (pregnancy) and changes in size (weight loss or gain, oedema) **require adaptation of solution size.**

Good to know

The quality of the socks' manufacture and finish is guaranteed for **6 months from the date of first use.**

In the event of significant morphological changes, and at each renewal, your patients' measurements must be repeated to ensure proper treatment follow-up.

INDICATIONS	RAL	1 10 - 15 mmHg	1 18 - 21 mmHg	2 23 - 32 mmHg	3 34 - 46 mmHg
	ASQUAL	2 15 - 20 mmHg	3 20 - 36 mmHg		
Prevention of functional disorders linked to the prolonged standing		✓			
Heavy legs		✓	✓		
C0 Restless legs		✓	✓		
Evening oedema			✓		
C1 Telangiectasia or reticular veins (< 3 mm)		✓	✓		
Following sclerotherapy for telangiectasia or reticular veins (< 3 mm)			✓		
C2 Varicose veins (≥ 3 mm)			✓	✓	
Following sclerotherapy or surgery for varicose veins (≥ 3 mm)			✓	✓	
C3 Treatment of chronic venous oedema				✓	
C4 Trophic disorders (pigmentation, eczema, lipodermatosclerosis, atrophie blanche, corona phlebectatica)				✓	✓
C5 Healed venous ulcer				✓	✓
C6 Open ulcer				✓	✓
Prevention and/or treatment of venous thrombosis (lower limb) and post-thrombotic syndrome				✓	✓
Prevention of venous thrombosis during long journeys			✓		
Prevention of venous thrombosis in pregnant women or post-partum			✓		
Chronic venous disease associated with pregnancy or the post-partum period				✓	✓
To obtain higher pressure by layering several compression devices		✓	✓		
Orthostatic hypotension			✓	✓	
Treatment of lymphoedema (maintenance phase)			✓	✓	✓


VENOFLEX *discover our adapted solutions*

Venoflex solutions meet your patients' needs with a variety of solutions adapted to everyday situations.

TAKING CARE OF YOUR APPEARANCE:

Collection

Textured, patterned products for those who love style.



- ♀ VENOFLEX **Fantaisie** · OGEE VENOFLEX **Incognito** ABSOLU · SWISS DOTS
- VENOFLEX **Kokoon** ABSOLU · LOZENGES
- ♂ VENOFLEX **City Confort** FIL D'ÉCOSSE

Discret

Discreet, elegant solutions that enhance your legs like a city product. Available in both transparent and semi-transparent.



- ♀ VENOFLEX **Incognito** ABSOLU VENOFLEX **Kokoon** ABSOLU
- VENOFLEX **Micro** VENOFLEX **Elegance**
- ♂ VENOFLEX **Elegance**

COMFORT FIRST:

Sensitive

Products with vegetal or technical fibers that provide softness and respect the skin.



- ♀ VENOFLEX **Fast** LIN VENOFLEX **Fast** COTON VENOFLEX **Fast** LAINE
- VENOFLEX **Simply** COTON FIN
- ♂ VENOFLEX **Fast** LIN VENOFLEX **Fast** COTON VENOFLEX **Fast** LAINE VENOFLEX **Fast** AIR
- VENOFLEX **City Confort** COTON VENOFLEX **Men**

SPECIFIC SITUATIONS:

Maternity

Tights specifically designed and adapted for mothers-to-be.



VENOFLEX **Kokoon** ABSOLU · MATERNITY VENOFLEX **Secret** VENOFLEX **Micro** VENOFLEX **Elegance**

Specific

Products specifically designed for medical situations and adapted to sensitive skin.



♀ VENOFLEX **Soft & Care** VENOFLEX **Clinic**

ACCESSIBLE SOLUTIONS*:

Original

A range that is simple.



- ♀ VENOFLEX **Secret** VENOFLEX **Secret** · MATERNITY VENOFLEX **Secret** OPAQUE
- ♂ VENOFLEX **Elégance**

* An entry-level range.

WOMEN

ASQUAL

**TRAVEL**

Socks
Easy to put on*
Comfort**

**LEISURE
AND WELLNESS**

Sensory comfort**
Sedentary or simple activities
Natural fibers

**ACTIVITIES** ++

Standing
Outdoor activities
Resistance

**CITY**

Professional activities
Evenings
Discretion

**PREGNANCY**

Special tights
for pregnant women

**JOINT
DISCOMFORT*****

Easy to put on
Arthritic hands


**PROLONGED
BED REST**

High-risk medical
and surgical situations

**STANDARD***

Simple solutions

Sensitive

VENOFLEX
Fast LINVENOFLEX
Fast COTONVENOFLEX
Fast LAINEVENOFLEX
Fast LINVENOFLEX
Fast COTONVENOFLEX
Fast LAINEVENOFLEX
Simply COTON FINVENOFLEX
Fast LINVENOFLEX
Fast COTONVENOFLEX
Fast LAINEVENOFLEX
Simply COTON FIN

Collection

VENOFLEX
Incognito ABSOLU • SWISS DOTSVENOFLEX
Fantaisie • OGEEVENOFLEX
Kokoon ABSOLU • LOZENGES

Discreet

VENOFLEX
Incognito ABSOLUVENOFLEX
Kokoon ABSOLU

Maternity

VENOFLEX
Kokoon ABSOLU • MATERNITY

Specific

VENOFLEX
Soft & Care

Specific

VENOFLEX
Clinic

Original

VENOFLEX
SecretVENOFLEX
Secret OPAQUEVENOFLEX
Secret • MATERNITY

* An entry-level range.

WOMEN

RAL



CITY

Professional activities
Evenings
Discretion

Discreet



VENOFLEX
Micro



VENOFLEX
Elegance

Maternity



VENOFLEX
Micro



VENOFLEX
Elegance



* FAST technology: easy to put on and take off.

** See product composition and properties.

*** Study EH2798 Tecni'sens for Thuasne®- February 2019. Ease of donning and doffing, comfort, softness and purchase intention of Thuasne® Soft & Care class 3 vs Sigvaris Expert, Varisma Comfort coton, Legger adiante 93 Coton and Varisan Elle & Lui, on 66 patients. Overall rating of 7.66/10 after use, significantly higher than other products tested. Intention to buy 71%.

MEN

ASQUAL

 TRAVEL Socks Easy to put on* Comfort**		 VENOFLEX Fast AIR
		 VENOFLEX Fast LIN
		 VENOFLEX Fast COTON
		 VENOFLEX Fast LAINE
 LEISURE AND WELLNESS Sensory comfort** Sedentary or simple activities Natural fibers		 VENOFLEX Fast AIR
		 VENOFLEX Fast LIN
		 VENOFLEX Fast COTON
		 VENOFLEX Fast LAINE
		   VENOFLEX City Confort COTON
 ACTIVITIES ++ Standing Outdoor activities Resistance		 VENOFLEX Fast LIN
		 VENOFLEX Fast COTON
		 VENOFLEX Fast LAINE
		   VENOFLEX City Confort COTON
 CITY Professional activities Evenings Discretion		   VENOFLEX City Confort COTON
	Collection	  VENOFLEX City Confort FIL D'ÉCOSSE
 JOINT DISCOMFORT*** Easy to put on Arthritic hands	Specific	 VENOFLEX Soft & Care
 PROLONGED BED REST High-risk medical and surgical situations	Specific	  VENOFLEX Clinic
 STANDARD* Simple solutions	Original	  VENOFLEX Elégance

* An entry-level range.

MEN

RAL



CITY

Professional activities
Evenings
Discretion

Discreet



VENOFLEX
Elegance

Sensitive



VENOFLEX
Men

VENOFLEX
Men



* FAST technology: easy to put on and take off.

** See product composition and properties.

*** Study EH2798 Tecni'sens for Thuasne® - February 2019. Ease of donning and doffing, comfort, softness and purchase intention of Thuasne® Soft & Care class 3 vs Sigvaris Expert, Varisma Comfort coton, Legger, Radiante 93 Coton and Varisan Elle & Lui, on 66 patients. Overall rating of 7.66/10 after use, significantly higher than other products tested. Intention to buy 71%.

Collection

WOMEN



Product worn: Venoflex Incognito Absolu Lozenges



City



Evening

CLASS

2

15 – 20 mmHg


Collection
WOMEN

Textured, patterned products for those who love style.

Wearing effective compression can also be about feminine style.
A range for trend-conscious women who attach great importance to their outfits.




Semi-transparent



Geometric & retro

VENOFLEX
Fantaisie • OGEE



Transparent



Light & timeless

VENOFLEX
Incognito **ABSOLU** • SWISS DOTS


Semi-transparent



Comfortable* & modern

VENOFLEX
Kokoon **ABSOLU** • LOZENGES

CLASS

2

VENOFLEX
Fantaisie • OGEE / p. 26

CLASS

2

VENOFLEX
Incognito **ABSOLU** • SWISS DOTS / p. 28

CLASS

2

VENOFLEX
Kokoon **ABSOLU** • LOZENGES / p. 30

* Use of a yarn that provides high elasticity and a ribbed cuff that does not create excess pressure on the upper calf.

These products are regulated Class I medical devices (MD) CE marked. Read instructions before use.

/ 25

For a retro geometric look

STOCKINGS

Breathable anti-slip**
for all-day wear.

TIGHTS

Stretchy panty
that won't compress*.

SUITABLE FOR:

- working days
- evenings
- mid-season

EASE OF DONNING

Elasticity of yarns used
in product construction.

RESISTANCE

Reinforced toe and sole
areas on all models.

COMPOSITION

Polyamide, elastane.

The stockings contain silicone.

Machine washable at  40°



VENOFLEX *Fantaisie* • OGEE

Feminine, retro compression



TRENDY

An ultra-feminine geometric pattern.



SEMI-TRANSPARENT

Medium mesh density combined with smooth, shiny yarns for a semi-transparent finish that highlights the legs.



GEOMETRIC & RETRO



www.asqual.com



STOCKINGS



TIGHTS

Reference

6325 01

6323 01




black
ogee

2

2

TIGHTS
Stretchy shaping panty that won't compress*.



STOCKINGS
Breathable anti-slip** for all-day wear.



RESISTANCE
Reinforced toe and sole areas on all models.

* Use of a thread for high elasticity.

** Zone constructed with 5 independent silicone waves.

Breathability, comfort** and femininity*

SOCKS

Patented bi-zone cuff with no compression:

- Single-layer zone with anti-slip thread: improves fit and breathability*.
- Double-layer zone: improved comfort.

STOCKINGS

Anti-slip wave effect

- Comfortable to wear.
- Good fit on the leg*: adapts to thigh morphology.

TIGHTS

Stretchy shaping panty with no pressure on the tummy, adapts to the shape of the body**.

SUITABLE FOR:

- working days
- evenings
- warmer seasons

EASE OF DONNING

A product that's easy to put on and take off*.
Stretchable, compression-free toes.

RESISTANCE

Unglazed toe. Reinforced zones at the ball of the foot and between the mesh and the anti-slip band*.

COMPOSITION

Polyamide, elastane.

The stockings contain silicone.

Machine washable at  40°C



VENOFLEX *Incognito* **ABSOLU** • SWISS DOTS

Feminine and timeless compression



TRENDY

Swiss dots,
the elegant,
timeless pattern.



TRANSPARENT

Low mesh density combined
with smooth, shiny yarns
for a sheer finish that sublimates
the legs like a city product.



LIGHT & TIMELESS



	Reference	
 SOCKS	5020 02	 black swiss dots
 STOCKINGS	5025 02	
 TIGHTS	5023 02	

TIGHTS
Stretchy shaping
panty with no pressure
on the tummy, adapts
to the shape of the
body**.



SOCKS

Patented bi-zone cuff with no compression.
- Single layer zone with anti-slip thread:
improves fit and breathability*.
- Double layer zone for improved **comfort**.



STOCKINGS

Anti-slip wave effect. **Comfortable
to wear and good fit on the leg***:
adapts to thigh morphology.



RESISTANCE

Reinforced toe and sole
areas on all models.

* Randomised, comparative, crossover clinical study, blind test with 32 women, reference DermScan 152E2280, January 2016

** Use of a yarn providing high elasticity.

Comfortable and assertive look

SOCKS

Bi-zone lining, stretchy ribbing*
limits the tourniquet effect**.

STOCKINGS

Breathable, non-slip silicone
"wave effect"***.

TIGHTS

Highly stretchable shaping panty
with no pressure on the tummy,
adapts to body shape*.

SUITABLE FOR:

- working days
- evenings
- colder seasons

EASE OF DONNING

Elasticity of yarns used
in product construction.

RESISTANCE

Reinforced toe and sole areas
on all models.

COMPOSITION

Polyamide, elastane.

The stockings contain silicone.

Machine washable at  40



VENOFLEX Kokoon

ABSOLU • LOZENGES



2
15 - 20 mmHg

Feminine and sophisticated compression



TRENDY

Diamond, the sophisticated pattern to assert your style.



SEMI-TRANSPARENT

Medium mesh density combined with smooth, shiny yarns for a semi-transparent finish that enhances the legs.



COMFORTABLE* & MODERN



	Reference	black diamonds
 SOCKS	5220 03	2
 STOCKINGS	5225 03	2
 TIGHTS	5223 03	2

TIGHTS

Highly **stretchable** shaping panty with no pressure on the tummy and adapts to body shape*.




SOCKS

Bi-zone lining, stretchy ribbing* limits the tourniquet effect**.



STOCKINGS

Breathable, non-slip silicone "wave effect"***.



RESISTANCE

Reinforced toe and sole areas on all models.



* The use of a yarn that is highly elastic.

** Ribbing that doesn't create excess pressure on the upper calf.

*** Focus group conducted in France and in Poland with 29 patients - Sensory test - OPINION WAY - October 2020.

Collection

MEN



Product worn: Venoflex City Confort Fil d'Écosse in black



Work



Elegance

CLASS

2

15 - 20 mmHg

Collection
MEN

Elegant, textured solutions for those who like to look stylish.

For everyday wear, solutions for active men looking for original products. Solutions designed for those who like to look their best and adapt their compression to their look.

CLASS

2

VENOFLEX

City Confort

FIL D'ÉCOSSE

/ p. 34

*Comfortable,
fine
and refined*

SUITABLE FOR:

- working days
- evenings
- warmer seasons

COMPOSITION

Cotton, polyamide, elastane.
The stockings contain silicone.

Machine washable at  40



VENOFLEX City Comfort

FIL D'ÉCOSSE

The refined compression



FINE AND ELEGANT

Fine, silky knit thanks to cotton yarn from Scotland*.



COMFORTABLE

Anatomically shaped to fit the calf and reduce the tourniquet effect.





	Reference	beige dakar	black	navy blue
 SOCKS	5329 02	2	2	2
 SOCKS LONG FOOT	5329 02	-	2	-
 SOCKS CALF +/-	5329 02	-	2	-
 OPEN-TOE SOCKS	5329 02	-	2	-
 STOCKINGS	5328 02	2	2	2
 OPEN-TOE STOCKINGS	5328 02	-	2	-



KNIT

Fine ribbed texture and silky feel thanks to the cotton yarn from Scotland*.



COMFORTABLE

An anatomical shape that hugs the contour of the calf to reduce the tourniquet effect.

* https://www.la-federation.com/data/club-maille_58014/categorie/105/charte_qualite_fr_462d7.pdf.

These products are regulated Class I medical devices (MD) CE marked. Read instructions before use.

Discreet WOMEN



Product worn: Venoflex Incognito Absolu in dark brown,
Venoflex Incognito Absolu in natural beige



Work



Elegance

CLASS

②

15 – 20 mmHg

③

20 – 36 mmHg

① RAL

18 – 21 mmHg

② RAL

23 – 32 mmHg

Discreet
WOMEN

*Discreet, feminine solutions
that enhance your legs like
a city product. Available in both
transparent and semi-transparent.*

For everyday wear, solutions for active women who want
to take care of their health with discretion.

CLASS ② VENOFLEX
Incognito ABSOLU / p. 38

CLASS ②
③ VENOFLEX
Kokoon ABSOLU / p. 42

CLASS ①
② RAL VENOFLEX
Micro / p. 44


CLASS ② RAL VENOFLEX
Elegance / p. 46

Discreet
WOMEN



VENOFLEX *Incognito* ABSOLU

*Range of transparent compression stockings,
ideal for a natural effect.*



A shade for different skin tones.

The range expands to offer **an inclusive compression solution** that matches **different skin tones**.

A range of colors from light to dark for a natural or tanned effect, depending on the skin tone.

*Hold, comfort** and breathability**

SOCKS:

Patented bi-zone cuff
with no compression:

- Single-layer zone with anti-slip thread:
improves fit/breathability*.
- Double-layer zone: improved comfort.

STOCKINGS:

Anti-slip wave effect

- Comfortable to wear*.
- Good fit on the leg*:
adapts to thigh morphology.

TIGHTS:

- Stretchy shaping panty
with no pressure on tummy,
adapts to body shape**.

SUITABLE FOR:

- working days
- warmer seasons

EASE OF DONNING

A product that's easy
to put on and take off*.
Stretchy, compression-free toes.

SOLIDITY

Run-resistant toe caps.
Reinforced zones on the sole
of the foot and between the mesh
and the anti-slip band*.

COMPOSITION

Polyamide, elastane.
The stockings contain silicone.

Machine washable at 



VENOFLEX

Incognito ABSOLU*Feminine, transparent compression***FEMINITY**

A wide palette of nude shades for spring and summer for a natural effect.

**TRANSPARENCY**

Low mesh density combined with smooth, shiny yarns for a transparent finish that sublimates the legs like a city product.

**SKIN-TONE NUDE PALETTE**

	Reference	natural	nude	golden	bronzed	dark brown	black
SOCKS	5020 02	2	2	2	2	2	2
OPEN-TOE SOCKS	5020 02	2	-	-	2	-	2
STOCKINGS	5025 02	2	2	2	2	2	2
OPEN-TOE STOCKINGS	5025 02	2	-	-	2	-	2
TIGHTS	5023 02	2	-	2	2	-	2

* Randomised, comparative, crossover clinical study, blind test with 32 women, reference DermScan 152E2280, January 2016.

** Use of a yarn providing high elasticity.

These products are regulated Class I medical devices (MD) CE marked. Read instructions before use.

All-day comfort and fit

SOCKS

Bi-zone lining, stretchy ribbing* limits
tourniquet effect**.

STOCKINGS

Breathable, non-slip silicone
"wave effect"***.

TIGHTS

Highly stretchable shaping
panty with no pressure on the tummy,
adapts to body shape*.

SUITABLE FOR:

- working days
- cold seasons

EASE OF DONNING

Elasticity of yarns used
in product construction.

RESISTANCE

Reinforced zones on all models
in the toes and sole of the foot.

COMPOSITION

Polyamide, elastane.

The stockings contain silicone.

Machine washable at  40°



VENOFLEX *Kokoon* ABSOLU

Soft, feminine compression



CLASS

②

15 – 20 mmHg

③

20 – 36 mmHg

Discreet
WOMEN



SOFT****

Microfiber effect thanks to the use of specific weave yarns.



SEMI-TRANSPARENT

Medium mesh density combined with smooth, shiny yarns for a semi-transparent finish that enhances the legs.

 www.asqual.com	Reference	 navy  honey  natural beige  black			
 SOCKS	5220 03	②	②	②	②
	5230 03	-	-	③	③
 SOCKS CALF +/-	5220 03	-	-	②	②
 OPEN-TOE SOCKS	5220 03	-	②	-	-
	5230 03	-	-	③	③
 STOCKINGS	5225 03	②	②	②	②
	5235 03	-	-	③	③
 STOCKINGS +/-	5225 03	-	-	②	②
 OPEN-TOE STOCKINGS	5225 03	-	②	-	-
	5235 03	-	-	③	③
 TIGHTS	5223 03	②	-	②	②
	5233 03	-	-	③	③

* The use of a yarn that is highly elastic.

** Ribbing that doesn't create excess pressure on the upper calf.

*** Focus group conducted in France and in Poland with 29 patients - Sensory test - OPINION WAY - October 2020.

These products are regulated Class I medical devices (MD) CE marked. Read instructions before use.

Comfort and flexibility**

- Morphological adaptation of the leg
- Special cuff and anti-slip
- Skin-friendly*

SUITABLE FOR:

- working days
- all seasons

EASE OF DONNING

Elasticity of yarns used
in product construction.

SOLIDITY

Reinforced heel.
Flat anti-slip seams
that reduces the risk
of tearing.

COMPOSITION

Polyamide, elastane.
The stockings contain silicone.

Machine washable at 



VENOFLEX Micro

Daily opaque compression

1

18 – 21 mmHg

2

23 – 32 mmHg



SOFT

Product made from microfibers for great softness when worn.



OPAQUE

High mesh density with textured yarn for an opaque finish.



SOCKS



STOCKINGS



STOCKINGS +



TIGHTS

Reference



light



caramel



tan
beige



black

30021A

1

1

1

1

30021B

2

2

2

2

30023D

1

1

1

1

30023E

2

2

2

2

30023G

1

1

1

1

30023H

2

2

2

2

30025A

1

1

1

1

30025B

2

2

2

2

Legal manufacturer Thuasne® Deutschland



Thuasne Deutschland GmbH
Sydney Garden 9
D-30539 Hannover

* Product certified Oeko-Tex 100.

** Use of a yarn providing high elasticity.

These products are regulated Class I medical devices (MD) CE marked. Read instructions before use.

Everyday comfort

- High-quality yarns
- Comfortable cuff and anti-slip
- Skin-friendly*

SUITABLE FOR:

- working days
- evenings
- all seasons

EASE OF DONNING

Elasticity of yarns used
in product construction.

RESISTANCE

Thanks to its particularly sturdy
construction, this round-knit
product is extremely resistant.

COMPOSITION

Polyamide, elastane.

The stockings contain silicone.

Machine washable at  40



VENOFLEX *Elegance*

Feminine, opaque compression



FEMINITY

A wide palette of discreet colors to match every wardrobe.



OPAQUE

High mesh density with textured yarn for an opaque finish.



SOCKS



STOCKINGS



TIGHTS

Reference	light	caramel	tan beige	diamond	navy	black
30021B	2	2	2	2	2	2
30023E	2	2	2	2	2	2
30025B	2	2	2	2	2	2

Legal manufacturer Thuasne® Deutschland



Thuasne Deutschland GmbH
Sydney Garden 9
D-30539 Hannover

* Product certified Oeko-Tex 100.

These products are regulated Class I medical devices (MD) CE marked. Read instructions before use.

Discreet

MEN



Product worn: Venoflex Elegance in black



Work



Elegance

CLASS

2 RAL

23 - 32 mmHg

Discreet
MEN

***Discreet solution
that looks like
a city product both
soft and elegant.***

For everyday wear, solution for active men who want to take care of their health with discretion.

CLASS **2** RAL | VENOFLEX
Elegance / p. 50

Everyday comfort

SUITABLE FOR:

- working days
- evenings
- all seasons

EASE OF DONNING

Elasticity of yarns used
in product construction.

RESISTANCE

Thanks to its particularly sturdy
construction, this round-knit
product is extremely resistant.

COMPOSITION

Polyamide, elastane.

The stockings contain silicone.

Machine washable at  40°C



VENOFLEX *Elegance*

Daily elegant compression



SOBRIETY

A wide palette of elegant colors to match every outfit.



COMFORT**

- High-quality yarns
- Comfortable cuff and anti-slip
- Skin-friendly*



	Reference	light	caramel	tan beige	diamond	navy	black
 SOCKS	20021B	2	2	2	2	2	2
 STOCKINGS	20023E	2	2	2	2	2	2
 STOCKINGS -/+	20023H	2	2	2	2	2	2
 TIGHTS	20025B	2	2	2	2	2	2

Legal manufacturer Thuasne® Deutschland



Thuasne Deutschland GmbH
Sydney Garden 9
D-30539 Hannover

* Product certified Oeko-Tex 100.

** Use of a yarn providing high elasticity.

These products are regulated Class I medical devices (MD) CE marked. Read instructions before use.

Sensitive

WOMEN



Product worn Venoflex Simply Coton Fin in powder pink



Technical
fibers



Comfort

CLASS

②

15 – 20 mmHg

③

20 – 36 mmHg

Sensitive
WOMEN

*Products with vegetal or technical fibers that provide softness and respect the skin.**

Solutions designed for women seeking sensorial comfort**.
Products that are both seasonally adapted and skin-friendly.

CLASS ② | VENOFLEX
Fast LIN / p. 54

CLASS ②
③ | VENOFLEX
Fast COTON / p. 56

CLASS ②
③ | VENOFLEX
Fast LAINE / p. 58

CLASS ② | VENOFLEX
Simply COTON FIN / p. 60

* Based on technical fibers. SENSIL® Aquarius facilitates absorption and evacuation of perspiration. <http://www.nilit.com> : SENSIL® BREEZE: soft fiber with cooling effect / SENSIL® AQUARIUS: high-performance hydrophilic fiber.

** Use of a yarn providing high elasticity.

Sensitive
WOMEN

*Natural,
breathable
and soft* comfort*

BREATHABLE

Thermoregulating natural linen fiber.

SOFTNESS

Combination of linen fiber
and viscose**.

SUITABLE FOR:

- traveling
- warmer seasons

COMPOSITION

Viscose, polyamide, linen, elastane.

Machine washable at  40°

Product worn: Venoflex Fast Lin
in natural stripe



VENOFLEX *Fast* LIN

Effortless and Lightweight**** compression*



New 30% thicker yarn for extra heel durability***.



LIGHTNESS

A thin, breathable and thermoregulating sock**.



FAST TECHNOLOGY FOR COMFORT* AND EASE OF DONNING

- Stitching on the top of the foot and no cuffs.
- Stretchy socks in height and width.



Reference

6120 04



azure blue

2



black strips

2



natural strips

2

Linen



Linen is a fiber that guarantees a resource-efficient and sustainable form of agriculture, relying solely on rainwater and sunlight. Linen is an inspiring fiber, thanks in particular to its **thermoregulating properties, breathable in summer and insulating in winter****.

The innovative and unique combination of **flax and modal provides exceptional softness** both inside and out, along with great suppleness and lightness, thereby enhancing therapeutic compliance, especially during the **summer months****.



Masters of
FLAX FIBRE

The certification guarantees sustainable farming practices and promotes the unique expertise of European fiber producers.

Masters of FLAX FIBRE™ certifies every stage of the transformation process, from fiber to finished products, intended for all markets.

* FAST technology: easy to put on and take off.

** Natural linen fiber with thermoregulating properties. <http://news.europeanflax.com/fr/lin/benefices-lin/>

*** Trial labo n°2023-356 made in July 2023.

**** Everything you need to know about European flax | ALLIANCE. (s. d.). <https://allianceflaxlinenhemp.eu/fr/tout-savoir-lin-europeen>.

These products are regulated Class I medical devices (MD) CE marked. Read instructions before use.

Sensitive
WOMEN

The daily comfort of cotton

SUITABLE FOR:

- outdoor activities
- all seasons

COMFORT

Insulating and absorbent
thanks to the cotton fiber.

COMPOSITION

Cotton, polyamide, elastane.

Machine washable at  40

Product worn: Venoflex Fast Coton
in grey irish marl



VENOFLEX *Fast* COTON

Effortless and resistant** compression*

2

15 - 20 mmHg

3

20 - 36 mmHg



RESISTANT**

Made of cotton to guarantee the strength of the mesh and increase the sock's durability.



FAST TECHNOLOGY FOR COMFORT* AND EASE OF DONNING

- Stitching on the top of the foot and no cuffs.
- Stretchy socks in height and width.



SOCKS

Reference

6120 01

6130 01

grey irish
marl

2

black
irish

2

-

3

Cotton



Cotton is valued for its properties that provide **strength and durability**. Cotton is well-known for its excellent moisture absorption, easy care, and strong resistance to UV rays. It is the fiber favored by **sensitive skin****.

* FAST technology: easy to put on and take off.

** https://www.la-federation.com/data/chiffres-cles_3ec82/fiche/633/fibres_naturelles_36d9d.pdf

These products are regulated Class I medical devices (MD) CE marked. Read instructions before use.

*Natural comfort,
insulation and
softness**.*

SUITABLE FOR:

- outdoor activities
- cold and humid weather

SOFTNESS

Knitted with Air-jet yarn for soft, supple knitwear.

ELEGANCE

A textured, elegant pattern thanks to the Jacquard knitting.

COMPOSITION

Wool, polyamide, elastane.

Machine washable at 



VENOFLEX *Fast* LAINE

Effortless and insulating** compression*

2

15 – 20 mmHg

3

20 – 36 mmHg



INSULATING AND THERMO-REGULATING**

Made from Merino wool to protect the legs from the cold and wick away perspiration.



FAST TECHNOLOGY FOR COMFORT AND EASE OF DONNING

- Stitching on the top of the foot and no cuffs.
- Stretchy socks in height and width.



SOCKS

Reference

6120 02

6130 02

natural
twistsanthracite
twistsblack
twists

2

2

2

3

-

-

*Merino wool***



The Venoflex Fast Laine innovation is made with 39% **merino wool**, designed to withstand all temperatures. **Three times finer than traditional wool**, merino wool contains a large amount of air, up to 80% of its own weight. **This gives it excellent insulating properties, acting as a natural barrier against the cold** while reducing the sensation of itchiness**.

* FAST technology: easy to put on and take off.

** https://www.la-federation.com/data/chiffres-cles_3ec82/fiche/633/fibres_naturelles_36d9d.pdf

These products are regulated Class I medical devices (MD) CE marked. Read instructions before use.

*Femininity,
resistance
and everyday
comfort*

FEMININITY:

- fine ribbed pattern
- wide range of soft colors.

SUITABLE FOR:

- leisure or travel
- mid-season

COMPOSITION

Cotton, polyamide, elastane.

The stockings contain silicone.

Machine washable at  40



VENOFLEX *Simply* COTON FIN

Comfortable and resistant compression



RESISTANT*

Made of cotton to guarantee mesh strength and increase durability.



COMFORTABLE

Anatomical shape that hugs the calf without compressing.



COLOURFUL



	Reference	natural beige	grey violet	powder pink	cream	black
SOCKS	5320 03	2	2	2	2	2
STOCKINGS	5325 03	2	-	-	-	2
TIGHTS	5323 03	-	-	-	-	2

Cotton



Cotton is valued for its properties that provide **strength and durability**. Cotton is well-known for its excellent moisture absorption, easy care, and strong resistance to UV rays. It is the fiber favored by **sensitive skin****.

* https://www.la-federation.com/data/chiffres-cles_3ec82/fiche/633/fibres_naturelles_36d9d.pdf

** An anatomical shape that hugs the shape of the calf without compressing.

These products are regulated Class I medical devices (MD) CE marked. Read instructions before use.

Sensitive

MEN



Product worn:
Venoflex Fast Laine in black



CLASS

1

10 – 15 mmHg

2

15 – 20 mmHg

3

20 – 36 mmHg

2 ^{RAL}

23 – 32 mmHg

Sensitive
MEN

Products with vegetal or technical fibers that provide softness and respect the skin.*

Solutions designed for men seeking optimal sensory comfort**.
Products that are both seasonally adapted and skin-friendly*.

CLASS 2 | VENOFLEX
3 | **Fast AIR** / p. 64

CLASS 2 | VENOFLEX
3 | **Fast LIN** / p. 66

CLASS 2 | VENOFLEX
3 | **Fast COTON** / p. 68

CLASS 2 | VENOFLEX
3 | **Fast LAINE** / p. 70

CLASS 1 | VENOFLEX
2 | **City Confort COTON** / p. 72
3

CLASS 2 | VENOFLEX
RAL | **Men** / p. 74

* Based on technical fibers. SENSIL® Aquarius facilitates absorption and evacuation of perspiration. <http://www.nilit.com> : SENSIL® BREEZE: soft fiber with cooling effect / SENSIL® AQUARIUS: high-performance hydrophilic fiber.

** Use of a yarn providing high elasticity.

Sensitive
MEN

*Breathable,
supple and
thermoregulating
comfort***

A unique combination of fibers
to adapt to all life's situations:
intensive work day, stamping,
traveling, physical activities...

SUITABLE FOR:

- outdoor activities
- traveling
- warmer seasons

COMPOSITION

Polyamide, elastane.

Machine washable at 

Product worn: Venoflex Fast Air in anthracite



VENOFLEX *Fast* AIR

*Breathable** and effortless* compression*

②

15 – 20 mmHg

③

20 – 36 mmHg



ABSORBENT

Nilit Aquarius* fiber facilitates absorption and evacuation of perspiration.



REFRESHING

Nilit Breeze* fiber provides a cooling effect.



FAST TECHNOLOGY: COMFORT AND EASE OF DONNING

- Stitching on the top of the foot and no cuffs.
- Stretchy socks in height and width.



SOCKS

Reference

6129 03

6139 03



navy



anthracite



black

②

②

②

-

-

③

Technical fibers



Venoflex Fast Air is a unique combination of two advanced technical fibers:

- **Sensil® Breeze, a soft and cooling fiber** that reduces sensations of heat by providing enhanced ventilation*** and flexibility.
- **Sensil® Aquarius, a high-performance hydrophilic fiber**** with an original triple-T micro-channel structure that **facilitates the absorption and evaporation of perspiration.**

* FAST technology: easy to put on and take off.

** <http://www.nilit.com>: SENSIL BREEZE: soft fiber with cooling effect / SENSIL AQUARIUS: high-performance hydrophilic fiber.

*** Sensil BREEZE fiber was externally tested and approved for its ability to lower body temperature by the textile research center in Belgium, Centexbel.

Sensitive
MEN

*Natural,
breathable
and soft comfort***

BREATHABLE

Thermoregulating
natural linen fiber.

SOFTNESS

Combination of linen
fiber and viscose**.

SUITABLE FOR:

- traveling
- warmer seasons

COMPOSITION

Polyamide, linen, elastane.

Machine washable at  40

Product worn:
Venoflex Fast Lin in grège



VENOFLEX *Fast* LIN

*Light** and effortless* compression*



New 30% thicker yarn for extra heel durability***.



LIGHTNESS

A thin, breathable and thermoregulating sock**.



FAST TECHNOLOGY FOR COMFORT* AND EASE OF DONNING

- Stitching on the top of the foot and no cuffs.
- Stretchy socks in height and width.



SOCKS

Reference

6129 04



azure blue

2



greige

2



natural

2



black

2

Linen



Linen is a fiber that guarantees a resource-efficient and sustainable form of agriculture, relying solely on rainwater and sunlight. Linen is an inspiring fiber, thanks in particular to its **thermoregulating properties, breathable in summer and insulating in winter****.

The innovative and unique combination of **flax and modal provides exceptional softness** both inside and out, along with great suppleness and lightness, thereby enhancing therapeutic compliance, especially during the **summer months****.



Masters of
FLAX FIBRE

The certification guarantees sustainable farming practices and promotes the unique expertise of European fiber producers.

Masters of FLAX FIBRE™ certifies every stage of the transformation process, from fiber to finished products, intended for all markets.

* FAST technology: easy to put on and take off.

** Natural linen fiber with thermoregulating properties. <http://news.europeanflax.com/fr/lin/benefices-lin/>

*** Trial labo n°2023-356 made in July 2023.

These products are regulated Class I medical devices (MD) CE marked. Read instructions before use.

Sensitive
MEN

The everyday comfort of cotton

SUITABLE FOR:

- outdoor activities
- all seasons

COMFORT

Insulating and absorbent
thanks to the cotton fiber.

COMPOSITION

Cotton, polyamide, elastane.

Machine washable at  40

Product worn:
Venoflex Fast Cotton in slate ribbed



VENOFLEX *Fast* COTON

Effortless and resistant** compression*

2

15 - 20 mmHg

3

20 - 36 mmHg



RESISTANT**

Made of cotton to guarantee the strength of the mesh and increase the sock's durability.



FAST TECHNOLOGY : COMFORT* AND EASE OF DONNING

- Stitching on the top of the foot and no cuffs.
- Stretchy socks in height and width.



SOCKS

Reference

6129 01

6139 01

slate
ribbed

2

-

black
ribbed

2

3

Cotton



Cotton is valued for its properties that provide **strength and durability**. Cotton is well-known for its excellent moisture absorption, easy care, and strong resistance to UV rays. It is the fiber favored by **sensitive skin****.

* FAST technology: easy to put on and take off.

** https://www.la-federation.com/data/chiffres-cles_3ec82/fiche/633/fibres_naturelles_36d9d.pdf

These products are regulated Class I medical devices (MD) CE marked. Read instructions before use.

Natural, insulating, soft comfort.*

Knitted with a special yarn*** for soft, supple knitwear.

SUITABLE FOR:

- outdoor activities
- cold and humid weather

SOFTNESS

Knitted with Air-jet yarn for soft, supple knits.

ELEGANCE

A textured, elegant pattern thanks to the Jacquard knitting.

COMPOSITION

Wool, polyamide, elastane.

Machine washable at  40°

Product worn:
Venoflex Fast Line in black



VENOFLEX *Fast* LAINE

*Insulating** and effortless* compression*

2

15 – 20 mmHg

3

20 – 36 mmHg



INSULATING AND THERMO-REGULATING**

Made from Merino wool to protect the legs from the cold and wick away perspiration.



FAST TECHNOLOGY: COMFORT* AND EASE OF DONNING

- Stitching on the top of the foot and no cuffs.
- Stretchy socks in height and width.



SOCKS

Reference

6129 02

6139 02



anthracite



black

2

2

-

3

Merino wool**



The Venoflex Fast Laine innovation is made with 39% **merino wool**, designed to withstand all temperatures. **Three times finer than traditional wool**, merino wool contains a large amount of air, up to 80% of its own weight. **This gives it excellent insulating properties, acting as a natural barrier against the cold** while reducing the sensation of itchiness**.

* FAST technology: easy to put on and take off.

** https://www.la-federation.com/data/chiffres-cles_3ec82/fiche/633/fibres_naturelles_36d9d.pdf

*** Knitted with Air-jet yarn

Sensitive
MEN

Everyday strength and comfort

SUITABLE FOR:

- leisure or traveling
- mid-season

COMPOSITION

Polyamide, elastane, cotton.

The stockings contain silicone.

Machine washable at  40°C



VENOFLEX City Comfort

COTON

CLASS

1

10 – 15 mmHg

2

15 – 20 mmHg

3

20 – 36 mmHg

The comfortable and resistant compression



RESISTANT

Made of cotton to guarantee the strength of the mesh and increase the sock's durability*.



COMFORTABLE

Anatomically shaped to fit the calf and reduce the tourniquet effect.



	Reference	granite	black	navy blue
SOCKS	5319 01	-	1	-
	5329 01	2	2	2
	5339 01	-	3	-
OPEN-TOE SOCKS	5339 01	-	3	-
STOCKINGS	5328 01	2	2	2
	5338 01	-	3	-
OPEN-TOE STOCKINGS	5338 01	-	3	-
TIGHTS	5323 01	-	2	-

Cotton



Cotton is valued for its properties that provide **strength and durability**. Cotton is well-known for its excellent moisture absorption, easy care, and strong resistance to UV rays. It is the fiber favored by **sensitive skin****.

* https://www.la-federation.com/data/chiffres-cles_3ec82/fiche/633/fibres_naturelles_36d9d.pdf

** An anatomical shape that hugs the shape of the calf without compressing.

These products are regulated Class I medical devices (MD) CE marked. Read instructions before use.

The everyday comfort of cotton

SUITABLE FOR:

- outdoor activities
- working days
- all seasons

COMFORT

Insulating and absorbent
thanks to the cotton fiber.

RESISTANT*

Made of cotton to guarantee
the strength of the mesh
and increase the sock's durability.

COMPOSITION

Cotton, polyamide, elastane.

Machine washable at  40



VENOFLEX Men

The invisible compression



ELEGANT

The sock is specifically designed to look like a classic city sock, discreet and elegant.



COMFORT** AND EASE OF DONNING

- Stitching on the top of the foot and no cuffs.
- Stretchy socks in height and width.



SOCKS

Reference

40021B



navy



black



* https://www.la-federation.com/data/chiffres-cles_3ec82/fiche/633/fibres_naturelles_36d9d.pdf

** Use of a yarn providing high elasticity.

These products are regulated Class I medical devices (MD) CE marked. Read instructions before use.

Maternity



Product worn: Venoflex Kokoon Absolu Maternity in black



CLASS

2

15 – 20 mmHg

1

RAL

18 – 21 mmHg

2

RAL

23 – 32 mmHg

Tights specifically designed and adapted for mothers-to-be.

A range of products for mothers-to-be who are looking for relief for their legs without compressing their tummies. A solution adapted to their changing bodies.

CLASS 2 | VENOFLEX
Kokoon ABSOLU • MATERNITY / p. 78

CLASS 1
2 | VENOFLEX
Micro / p. 80

CLASS 2 | VENOFLEX
Elegance / p. 82

Also find

Original
MATERNITY

The Secret range, with its simplest maternity solution, is available in the "Original" segment.

See p. 97.

Comfort and hold all day long*

Shaping panties with high stretch
and no pressure on the tummy,
adapting to the body shape*.

SUITABLE FOR:

- rest
- everyday activities
- cold seasons

EASE OF DONNING

Stretchy non-slip
and specific ankle construction.

RESISTANCE

Reinforced toe and sole zones
on all models.

COMPOSITION

Polyamide, elastane.

Machine washable at 

Product worn: Venoflex Kokoon Absolu Maternity
in black diamonds



VENOFLEX Kokoon

ABSOLU • MATERNITY

*Soft**, feminine compression
designed for mothers-to-be.*



COMFORTABLE AND PRACTICAL

Compression-free, tummy-hugging stretchy panties in a semi-opaque finish that takes care of pregnant women without compromising elegance.



SOFT**

Microfiber effect thanks to the use of special yarns, designed for the mother and her baby.



MATERNITY
TIGHTS

Reference

5224 03



black
diamonds

2



natural
beige

2




black

2

TIGHTS

Stretchable,
pressure-free,
that gently follows the
shape of the belly and
adapts to the body's
morphology*.



RESISTANCE

Reinforced toe and sole areas
on all models.

* Use of a yarn that is highly elastic.

** Focus group conducted in France and in Poland with 29 patients - Sensory test - OPINION WAY - October 2020.

These products are regulated Class I medical devices (MD) CE marked. Read instructions before use.



Comfort** and flexibility

- Morphological adaptation of the leg
- Special cuff and anti-slip
- Skin-friendly*

SUITABLE FOR:

- working days
- all seasons

EASE OF DONNING

Elasticity of yarns used
in product construction.

SOLIDITY

Reinforced heel. Flat anti-slip seams
that reduce the risk of tearing.

COMPOSITION

Polyamide, elastane.

Machine washable at 

VENOFLEX Micro

Daily opaque compression

1

18 – 21 mmHg

2

23 – 32 mmHg



COMFORTABLE AND PRACTICAL

Compression-free, tummy-hugging stretchy panties in an opaque finish that takes care of pregnant women without compromising elegance.



SOFT

Product made from microfibers for great softness when worn.



MATERNITY
TIGHTS

Reference

30027A

30027B



light



caramel

tan
beige

black

1

1

1

1

2

2

2

2

Legal manufacturer Thuasne® Deutschland



Thuasne Deutschland GmbH
Sydney Garden 9
D-30539 Hannover

* Product certified Oeko-Tex 100.

** Use of a yarn providing high elasticity.

These products are regulated Class I medical devices (MD) CE marked. Read instructions before use.

Everyday comfort

- High-quality yarns
- Comfortable cuff and anti-slip
- Skin-friendly*

SUITABLE FOR:

- working days
- evenings
- all seasons

EASE OF DONNING

Elasticity of yarns used in product construction.

RESISTANCE

Thanks to its particularly sturdy construction, this round-knit product is extremely resistant.

COMPOSITION

Polyamide, elastane.

Machine washable at  40

Product worn: Venoflex Elegance in black



VENOFLEX *Elegance*

Feminine, opaque compression



FEMINITY

A wide palette of discreet colors to match every wardrobe.



COMFORTABLE AND PRACTICAL

Compression-free, tummy-hugging stretchy panties in an opaque finish that takes care of pregnant women without compromising elegance.



MATERNITY
TIGHTS

Reference

20027B



light



caramel



tan
beige



diamond



navy



black

2

2

2

2

2

2



Specific

UNISEX





Hospital



Comfort

CLASS

1

10 – 15 mmHg

2

15 – 20 mmHg

3

20 – 36 mmHg

Specific

Products specifically designed for medical situations and adapted to sensitive skin.*

Solutions conceived and designed for use in hospital environments, to respond to your patients' problems of joint discomfort and sensitive skin.

VENOFLEX *Soft & Care*



*Suitable for sensitive skin**

VENOFLEX *Clinic*



Prolonged bed rest

CLASS 2
3 | VENOFLEX
Soft & Care / p. 86

CLASS 1
2 | VENOFLEX
Clinic / p. 88

* https://www.la-federation.com/data/club-maille_58014/fiche/108/leaflet_fibres_synthetique_fr_080d5.pdf.

These products are regulated Class I medical devices (MD) CE marked. Read instructions before use.

***The Venoflex
solution,
best suited
to sensitive skin*.***

Viscose predominantly
present in socks*.


SUITABLE FOR:

- joint discomfort
- sensitive skin
- need for warmth

COMPOSITION

Viscose, polyamide, elastane.

Machine washable at  40



VENOFLEX Soft & Care

Compression for sensitive skin*.



CLASS

2

15 – 20 mmHg

3

20 – 36 mmHg

GENTLE FOR 94% OF
OF PATIENTS**



SOCKS FOR SENSITIVE SKIN*.

Made from viscose, a breathable plant-based fiber renowned for its softness**.



EASY TO PUT ON AND TAKE OFF FOR 91% OF PATIENTS**.

Thicker, easier-to-grip socks with expanded knitting.



SOCKS


Reference	
6220 01	2
6230 01	3



black

WARM AND COCOONING FEEL

Ideal for warming up the extremities, often cold due to a sedentary lifestyle*.



KNIT AND CUFF

- Made with **viscose** (a naturally derived fiber) for softness** and comfort for **sensitive skin**.
- **Wider knit: easier donning**, suitable for people with **osteoarthritis***.

* https://www.la-federation.com/data/club-maille_58014/fiche/108/leaflet_fibres_synthetique_fr_080d5.pdf

** Techni'sens EH2798 study for Thuasne® - February 2019. Ease of fitting and removal, comfort, softness and purchase intention for the Thuasne® Soft & Care class 3 product vs Sigvaris Expert, Varisma Confort coton, Legger, Radiante 93 Coton and Varisan Elle & Lui, on 66 patients. Our overall rating of 7,66/10 after use, is significantly higher than other products tested. 71% purchase intention.

These products are regulated Class I medical devices (MD) CE marked. Read instructions before use.

Compression for prolonged bed rest

SUITABLE FOR:

- bed rest
- high-risk operations.

PRECAUTIONS

When using an electric scalpel, it is advised to use a sock rather than a stocking to keep the device away from the scalpel. The sock should be far enough away so as not to make contact with the scalpel plate or electrodes.

COMPOSITION

Polyamide, elastane.

The stockings contain silicone.

Machine washable at  40°C



VENOFLEX Clinic

Compression for prolonged bed rest

1

10 – 15 mmHg

2

15 – 20 mmHg



HOSPITAL

Range for medical and surgical situations at risk of venous thrombosis (infection, prolonged bed rest, cancer, surgery...) or for pre-operative periods.



ANTISTATIC

These products do not contain metal wires that could interfere with a magnetic field.



www.asqual.com

Reference



white



ANTI-THROMBOSE
SOCKS

5820 01

2



ANTI-THROMBOSE
STOCKINGS

5815 01



1

5825 01

2

THIGH-HIGH STOCKINGS

By wearing Venoflex Clinic stockings, people who must undergo **prolonged bed rest** protect themselves from the risk of oedema or venous thrombosis.



KNIT/ANTISTATIC

Does not contain **metallic threads** that could interfere with a magnetic field.

OPEN TOE

- helps improving the air circulation at the foot;
- reduces discomfort at the toes;
- allows easy monitoring of distal pressure.



Original WOMEN



Product worn: Venoflex Secret Opaque in black



Accessible

CLASS

1

10 – 15 mmHg

2

15 – 20 mmHg

3

20 – 36 mmHg

Original
WOMEN

A basic range

For those who prefer a simple, basic solution.

CLASS 1
2
3

VENOFLEX
Secret / p. 92

CLASS 2

VENOFLEX
Secret **OPAQUE** / p. 94

Original
WOMEN

*Original:
effective and simple
compression*

SUITABLE FOR:

- work and daily activities
- warmer seasons

COMPOSITION

Polyamide, elastane.

The stockings contain silicone.

Machine washable at  40°



VENOFLEX Secret

Simple compression for women



SIMPLE

The classic compression solution available in socks, stockings and tights.



SEMI-TRANSPARENT

Medium mesh density combined with flat, shiny yarns for a semi-transparent finish.



	Reference	NEW						
		navy	natural beige	dune	golden	tan beige	dark brown	black
SOCKS	5910 01	-	-	-	1	-	-	-
	5920 01	2	2	2	2	2	2	2
	5930 01	-	-	-	3	-	-	3
SOCKS CALF +/-	5920 01	-	-	-	2	-	-	2
OPEN-TOE SOCKS	5920 01	-	-	-	2	2	2	2
	5930 01	-	-	-	3	-	-	3
STOCKINGS	5915 02	-	-	-	1	-	-	-
	5925 04	2	2	2	2	2	2	2
	5935 01	-	-	-	3	-	-	3
STOCKINGS +/-	5925 04	-	-	-	2	-	-	2
OPEN-TOE STOCKINGS	5925 04	-	-	-	2	2	2	2
	5935 01	-	-	-	3	-	-	3
TIGHTS	5923 01	-	-	-	2	-	-	2



*New design approved by patients**



Appearance:

- ☒ Soft to the touch
- ☒ More transparent**



When worn:

- ☒ Easier to remove**
- ☒ Soft feel on the skin
- ☒ Resistant and comfortable

*Focus group conducted in France with 18 female patients – Kantar – December 2023.

**Versus the previous version of Venoflex Secret.

Original:
WOMEN

*Original:
effective and simple
compression*

SUITABLE FOR:

- work and daily activities
- cold seasons

COMPOSITION

Polyamide, elastane.

The stockings contain silicone.

Machine washable at  40°C

Product worn: Venoflex Secret Opaque in black



VENOFLEX *Secret* OPAQUE

Simple compression for women



SIMPLE

The classic compression solution available in socks, stockings and tights.



OPAQUE

High mesh density with textured yarn for an opaque finish.



	Reference	 navy	 golden	 black
 SOCKS	5920 05	2	2	2
 STOCKINGS	5925 05	2	2	2
 TIGHTS	5923 05	2	2	2

Original MATERNITY



Product worn: Venoflex Secret Maternity in golden



CLASS

2

15 - 20 mmHg

Original
MATERNITY

A basic range for mothers-to-be

For those who prefer a simple, basic solution.

CLASS 2

VENOFLEX

Secret • MATERNITY / p. 98

Also find

Maternity

A selection of maternity tights, with their premium solution, is also available in the "Maternity" segment.

See p. 77.

Original
MATERNITY

Original: effective and simple compression

SUITABLE FOR:

- rest
- everyday activities
- warm seasons

COMPOSITION

Polyamide, elastane.

Machine washable at  40°

Product worn: Venoflex Secret Maternity in golden



VENOFLEX *Secret* • MATERNITY

Simple compression for mothers-to-be



SIMPLE

The classic compression solution designed for pregnant women.



SEMI-TRANSPARENT

Medium mesh density combined with flat, shiny yarns for a semi-transparent finish.



MATERNITY
TIGHTS

Reference

5924 01



golden




black

2

2

TIGHTS
Stretchy
and pressure-free,
that hug the shape
of the belly and adapt
to the body's shape***.



RESISTANT

Reinforced toe and sole
areas on all models.



*New design approved by patients**



Appearance:

- ☒ Soft to the touch
- ☒ More transparent**



When worn:

- ☒ Easier to remove**
- ☒ Soft feel on the skin
- ☒ Resistant and comfortable

*Focus group conducted in France with 18 female patients - Kantar - December 2023.

**Versus the previous version of Venoflex Secret.

***Use of a yarn that is highly elastic.

Original MEN



Product worn: Venoflex Éléance in black



Accessible

CLASS

1

10 – 15 mmHg

2

15 – 20 mmHg

3

20 – 36 mmHg

Original
MEN

A basic range

For those who prefer a simple, basic solution.

CLASS
1
2
3

VENOFLEX
Élégance / p. 102

Original
MEN

Original: effective and simple compression

SUITABLE FOR:


- work and daily activities
- warm seasons

COMPOSITION

Polyamide, elastane.

The stockings contain silicone.

Machine washable at  40



VENOFLEX Élégance

Simple compression for men

CLASS

1

10 – 15 mmHg

2

15 – 20 mmHg

3

20 – 36 mmHg

Original
MEN



SIMPLE

Everyday compression.



EFFECTIVE


The classic compression solution available in socks and stockings.

 COMPRESSION MÉDICALE www.asqual.com		Reference	 beige sand	 navy	 grey	 black
 SOCKS	5119 01	-	-	-	1	
	5129 01	2	2	2	2	
	5139 01	3	-	-	3	
 SOCKS LONG FOOT	5129 01	-	2	-	-	
 SOCKS CALF +/-	5129 01	-	-	-	2	
 OPEN-TOE SOCKS	5139 01	3	-	-	3	
 STOCKINGS	5128 02	2	-	-	2	
	5138 01	3	-	-	3	
 OPEN-TOE STOCKINGS	5138 01	3	-	-	3	

Stocking puller

Ref. W2500001001

- Place the stocking puller on a table to make it easier to put on your sock or stocking.



1

- Present the stocking with the heel facing in your direction.




2

- Fold the upper edge of the stocking outward.




3

- Slide your stocking on the device until the heel appears.




4

- Put the stocking puller on the floor.
- Insert your foot inside.




5

- Take the stocking puller by the handles and gently pull up the stocking just below the knee with a back-and-forth movement.
- Push the stocking puller backward to put it aside.



6

- Finish applying the stocking manually without wrinkling or pulling.



7

- Massage your leg to finish smoothing out the stocking.

Slipper

Included in the ASQUAL class 3 products*

Ref. 23374508001

Aiding device to help putting on Thuasne® compression/retention products
(lower limbs)



CONVENIENT

2-in-1 device adapted for use with the open-toe AND closed-toe compression solutions.

LIGHT AND COMPACT

Packaged in its zipped pouch, patients can have it at hand in all circumstances.

WASHABLE

Machine washable at  40°C



*Except FAST and Soft & Care.

THUASNE® BANDAGES



Which prescription for which pathology?

	Biflex BANDAGES			MULTI TYPE KIT	SPECIFIC BANDAGES		
PRODUCT TYPE	BIFLEX 16+ Half-square overlap	BIFLEX 16+ Full-square overlap	BIFLEX 17+ Half-square overlap	KIT BIFLEX	COTTON BANDAGE	BIFLEXIDEAL < 0,9 et > 1.3	MOBIDERM® (Pad and bandages)
Application/use	Compression: immobile and mobile patients IPS above 0.9 and below 1.3				Skin protection sub-layer/strict retention of the limb in addition to compressive layer	Retention: mobile patients only	Drainage of oedemas under reduction bandage
CEAP CLASSIFICATION							
Varicose veins (C2): 15 - 36 mmHg	●					●	
Post-sclerotherapy/post-varicose surgery (C2): 15 - 36 mmHg	●					●	
Venous oedemas (C3): 20 - 36 mmHg	●			●	●	●	●
Pigmentation, venous eczema (C4a): 20 - 36 mmHg	●				●	●	
Venous hypodermatitis, atrophie blanche (C4b): 20 - 36 mmHg	●				●	●	
Corona phlebectatica (C4c): 20 - 36 mmHg	●				●	●	
Healed venous ulcer (C5): 20 - 36 mmHg or > 36 mmHg		●			●	●	
Active venous ulcer (C6):> 40 mmHg			●	●	●	●	
VENOUS THROMBOSIS: PREVENTION AND COMPLICATIONS 20 TO 36 mmHg							
Prevention of deep venous thrombosis and risk of relapse: 20 - 36 mmHg	●						
Superficial venous thrombosis lower limb: 20 - 36 mmHg	●					●	
Deep venous thrombosis lower limb: 20 - 36 mmHg	●					●	
Deep venous thrombosis upper limb: 20 - 36 mmHg	●					●	
Post-thrombotic syndrome: prevention as an alternative to compression stockings	●	●	●			●	
LYMPHOEDEMA							
Lymphoedema of the upper limb: use the maximum tolerated pressure	●	●	●	●	●	●	●
Lymphoedema of the lower limb: use the maximum tolerated pressure	●	●	●	●	●	●	●

HAS report/medical compression device for individual use: use in vascular diseases, september 2010

The pressure to consider is the highest tolerated by the patient.

Contra-indications

Biflex: PAD (peripheral arterial disease) with ABI < 0.6, skin wounds, abscess, infected skin conditions, extra-anatomic bypass, uncontrolled heart failure, phlegmasia cerulea dolens (painful blue oedema with arterial compression), septic thrombosis, hypersensitivity to any of the components, severe peripheral neuropathy.

Biflexideal: Do not use in the event of known allergy to any of the components, PAD (peripheral artery disease) with an ABI < 0.6, skin wounds, abscess, infected skin conditions of a limb concerned, advanced diabetic microangiopathy (if pressure > 30 mmHg), decompensated heart failure, phlegmasia cerulea dolens (painful blue phlebitis with arterial compression) of a limb concerned, extra-anatomic bypass, septic thrombosis of a limb concerned, severe peripheral neuropathy.

Kit Biflex: Arterial pathologies: peripheral arterial disease (PAD), arterial ulcer or mixed ulcer with a predominantly arterial component brachial index ABI < 0.6 (according to the French Health Authority) or > 1.3). Patients

suffering from advanced diabetic microangiopathy phlegmasia cerulea dolens with arterial compression, or septic thrombosis. Extra-anatomic bypass. Uncontrolled heart failure. Pre-existing or infected dermatological conditions (abscesses, boils, infected dermatitis). Known allergy to one of the components. Advanced peripheral neuropathy. Uncontrolled diabetes. Specifically when worn on an upper limb: brachial plexus pathology. Vasculitis of the extremities.

MOBIDERM®: Do not use on skin presenting infections of the limb or acute inflammation, known allergies to any of the components, septic thrombosis, severe peripheral neuropathy of the limb. Lower limbs: lower limb peripheral arterial disease (PAD) with ankle brachial index (ABI) < 0.6, decompensated heart failure, advanced diabetic microangiopathy, phlegmasia cerulea dolens (painful blue oedema with arterial compression), extra-anatomic bypass. Upper limbs: brachial plexus pathology, vasculitis of the extremities. The table above is only for reference when delivering a device. The level of pressure from the Biflex bandages varies according to the application technique (bandage tension, number of overlays...) and the patient's anatomic profile.



Kit Biflex

SINCE 1932

Elastic compression bandages

for the treatment of venous ulcers
of the lower limb.

(CEAP classification Stage C6)

- Chronic venous oedema
- Lymphoedema

WASHING PICTOGRAM

- Number of washes
- Washing temperature

x 6 SELF-GRIPPING FASTENERS


- Ensure a good maintenance of the strips
- Don't hurt the patient

x 1 FLEECY ELASTIC BANDAGE

- Protects the skin and bony prominences
- Contentive effect

x 1 TWO-WAYS ELASTIC BANDAGE

- Provides the level of pressure required on venous ulcer treatment, chronic venous oedema and lymphoedema



1 WIDTH: 10 CM

For an easier heel grip


AVAILABLE IN 3 SIZES AND 1 WIDTH TO FIT PATIENTS' MORPHOLOGIES:

Size	Ankle circumference	Bandage 1 length	Bandage 2 length	Bandages width	Reference
1	< 25 cm	4 m	3.5 m	10 cm	170050 20199999
2	25 - 30 cm	4.5 m	4 m		170050 20299999
3	> 30 cm	5 m	4.5 m		170050 20399999

HEALTH ALL SQUARED UP!

Kit Biflex

Plus points



REPOSITIONABLE

- The two bandages in the kit can be put back in place in the event of dressing changes or marked variations in oedema



WASHABLE

- The two bandages in the kit are machine washable at 40°C

Can be left in place day and night, for up to 7 days.

The 1st washable, removable and reusable Bi-layer Kit



2X MORE ULCERS HEALED AT 16 WEEKS⁽¹⁾

- Demonstrated Clinical Superiority



EASY TO USE⁽¹⁾

- Numbered and calibrated bandages with recovering lines



IMPROVEMENT OF THE PATIENT'S QUALITY OF LIFE⁽¹⁾

- Softness, breathability, comfort⁽²⁾
- Latex free

Washable and reusable bandages

-47% of Biflex Kits Used*

97.6% of patients are satisfied to very satisfied⁽¹⁾

97.8% of caregivers are satisfied to very satisfied⁽¹⁾

*Vs PROFORE® over a 16-weeks treatment period. Results on FAS population.

(1) Jean L. Gillet et al.: Clinical superiority of an innovative two-component compression system versus four-component compression system in treatment of active venous leg ulcers: A randomized trial; SAGE Publications 2019. KITBIFLEX is a class I medical device. It is indicated in the treatment of stage C6 venous ulcer, chronic venous oedema and lymphoedema. Read the leaflet carefully before use. Contraindications: Arterial pathologies, peripheral arterial disease (PAD), arterial ulcer or mixed ulcer with a predominantly arterial component (ankle-brachial index: ABI < 0.6 (according to the French Health Authority) or > 1.3; patients suffering from advanced diabetic microangiopathy, phlegmasia cerulea dolens, with arterial compression, or septic thrombosis, extra-anatomic bypass, uncontrolled heart failure, pre-existing or infected dermatological conditions (abscesses, boils, infected dermatitis), known allergy to one of the components, advanced peripheral neuropathy, uncontrolled diabetes.

(2) Compulse study: Evaluation of the clinical efficacy and safety of the KitBiflex 2-layer compression system versus the PROFORE® 4-layer compression system in the treatment of leg ulcers of venous origin (stage C6). National, multicenter, randomized, open-label, two parallel-group, phase IIIb clinical study. Study #: NCT02782689.

Compression bandages and accessories

Cutaneous protection



PROTECTIVE FOAM

Accessories for medical compression devices. Composition 100% cotton.

Width	Length	Color	Reference
11 cm	4 m	White	3111790230182

PROTECTIVE UNDERWRAP

Thin stretchy foam bandage

Protection of skin surface.

Width	Length	Color	Reference
7 cm	27 m	White	3111790258612

Containment

BIFLEXIDEAL

Short-stretch compression bandage

INDICATIONS: Treatment of oedema (traumatic or chronic) - Treatment of lymphoedema - Chronic venous insufficiency (as an alternative to compression stockings - Following sclerotherapy or varicose vein surgery - Treatment of venous ulcers (with SPI > 0.8 or under close medical supervision for SPI between 0.5 and 0.8) - Symptomatic treatment of acute, deep or superficial venous thrombosis - Prevention of post-thrombotic syndrome: Alternative to compression stockings for at-risk patients - Alternative to compression stockings for the treatment of leg oedema or swelling during pregnancy - Treatment of symptoms related to soft tissue injury. 65% extension at a force of 45 cN/cm.

Width	Length	Color	Reference
3 cm	5 m	White	3111790232216
6 cm			3111790232223
8 cm			3111790232230
10 cm			3111790232247

BIPLAST

Adhesive short-stretch compression bandage

INDICATIONS: Soft tissue or musculoskeletal system injuries. Strapping: post-trauma support. Fractures of the small bones (finger, foot, etc.). Congenital or acquired bone deformity of the foot. Vascular or lymphatic diseases requiring short-stretch or multicomponent bandage. Multi-directional adhesion to the edge of the tape. Flexible and adjustable. 50.2% maximum stretch

Width	Length	Color	Reference
3 cm	2,5 m	White	3111790072362
6 cm			3111790072379
8 cm			3111790072386
10 cm			3111790072393
15 cm			3111790072409
20 cm			3111790072416

Order in
France



SOFT ON THE SKIN⁽¹⁾

- designed to provide comfort
- can be used directly on the skin, providing protection under the bandage



TARGETED APPLICATION

- can be used in association to Thuasne® compressive bandages, Mobiderm® and Biflex Self Adjust



TARGETED APPLICATION AND THIN & PROTECTIVE

- can be used in association to Thuasne® compressive bandages, Mobiderm® and Biflex Self Adjust
- provides a thin layer of protection



SOFT ON THE SKIN⁽²⁾

- designed to provide comfort and breathability⁽³⁾



SHORT-STRETCH ELASTIC BANDAGE

- can be worn continuously, both day and night⁽⁴⁾
- provides good pressure retention over time⁽⁵⁾



WASHABLE AND REPOSITIONABLE⁽⁶⁾

- allows repositioning in case of a dressing change or significant variations in oedema
- washable 30 times⁽⁷⁾ at 40°



HERRINGBONE GLUING

- facilitates application and maintains flexibility of movement

(1) 100 % Cotton
(2) Soft on the skin: cotton + viscose (étude IFTH n°DEV004842)
(3) Designed to be comfortable: soft and breathable materials Flexibility: étude IFTH n° DEV004842.
(4) Can stay in place day and night (short stretch bandage) > low resting pressure: <https://journals.cambridge.org/ajw/volume-25-number-1/benefits-inelastic-and-short-stretch-bandages-leg-ulcer-narrative-review>
(5) Good pressure retention over time > elastic bandages : CER
(6) Reusable / repositionable: woven strips, no surface treatment (product characteristics).
(7) Stability test R&D.

Biflex

SINCE 1932

• Two-ways elastic bandage

Length elasticity/extensible width

• Two types of bandages:

BIFLEX 17+



BIFLEX 16+



WASHABLE AND REPOSITIONABLE⁽¹⁾

- can be put back in place in case of dressing changes, or significant variations in oedema
- washable 90 times⁽²⁾ at 40°



LONG-STRETCH ELASTIC BANDAGE

- allows good pressure maintenance over time, including during activity⁽³⁾
- adapted especially for inactive, hospitalized, bedridden or obese patient⁽⁴⁾



CALIBRATED⁽⁵⁾

- ensures the right level of pressure and therefore the effectiveness of the treatment
- provides a uniform and regular pressure level, appropriate to the limb
- allows the reproducibility of the application

Indications

- Alternative to compression stockings (15 - 36 mmHg) for varicose veins
- Chronic oedema
- Post-traumatic oedema
- After sclerotherapy or surgery of varicose veins as an alternative to compression stockings
- Alternative to compression stockings for trophic disorders (dermatitis/lipodermatosclerosis)
- Venous leg ulcers (multitype bandaging, elastic bandage for patients with reduced mobility)
- Symptomatic treatment of acute superficial and deep venous thrombosis
- Alternative to compression or antithrombotic stockings
- Post-thrombotic syndrome prevention: alternative to compression stockings for high-risk patients
- Lymphoedema treatment

CHOOSING ONE'S BIFLEX BANDAGE

Average height 40 cm

	Ankle	Below the knee	Top of thigh
	cB < 26 cm	cC < 39 cm	cG < 62 cm
Half-square overlap	8 cm x 1.5 m	8 cm x 3 m or 10 cm x 3 m	10 cm x 5 m
Full-square overlap	10 cm x 1.5 m	8 cm x 4 m or 10 cm x 3.5 m	8 cm x 5 m + 8 cm x 4 m or 10 cm x 5 m + 10 cm x 3 m
	cB ≥ 26 cm	cC ≥ 39 cm	cG ≥ 62 cm
Half-square overlap	10 cm x 1.5 m	8 cm x 4 m x 10 cm x 3.5 m	8 cm x 5 m + 8 cm x 3 m or 10 cm x 3.5 m + 10 cm x 3 m
Full-square overlap	8 cm x 3 m	8 cm x 5 m or 10 cm x 4 m	10 cm x 5 m + 10 cm x 4 m

(1) Reusable / repositionable: woven strips, no surface treatment (product characteristics).

(2) Stability test R&D.

(3) Allows good pressure maintenance over time, including during activity > elastic bandage.

(4) Can be worn on an inactive patient > patient profile: especially for hospitalized, bedridden, immobile or obese patients, alone or as a component of a multitype band.

(5) Advantages of calibration > patented calibration system: <https://zibocaredenmark.dk/wp-content/uploads/2021/09/UrgoK2-Hanna-study-British-Journal-of-Nursing.pdf> (analogie calibration Kit Biflex).

Compression bandages



Order in
France

Compression

BIFLEX+ GRADUATED

Graduated elastic compression bandage

Two-ways elastic bandage. Elastic in length and extensible in width. Repeatable control of the tension exerted thanks to the fitting indicator. Multiple overlappings possible to achieve multiple pressure level.



2 versions


	Color	Width	Length	Reference
Biflex 16⁺ (Light) 	Beige	8 cm	3 m	3111790234012
			3.5 m	3111790234029
			4 m	3111790234036
			5 m	3111790234043
		10 cm	3 m	3111790234067
			3.5 m	3111790234074
			4 m	3111790234081
			5 m	3111790234098
Biflex 17⁺ (Strong) 	Beige	8 cm	4 m	3111790234296
			3 m	3111790234302
		10 cm	3.5 m	3111790234319
			4 m	3111790234326
			5 m	3111790234333





BIFLEX APPLICATION IN 2 SIMPLE STEPS



1/ STRETCHING



2/ OVERLAPPING





Compression

BIFLEX+ PRATIC **Graduated elastic compression bandage with stirrup** Fitting stirrup makes this bandage easy to use. Enhances patient autonomy.



Designation	Color	Width	Length	Reference
Biflex16⁺Pratic (Light)	Beige	8 cm	3 m	160003 00830002
			4 m	160003 00840002
	Beige	10 cm	3 m	160003 01030002
			4 m	160003 01040002
Biflex17⁺Pratic (Strong)	Beige	10 cm	3 m	170003 01030002
			4 m	170003 01040002




PRATIC STIRRUP APPLICATION

1




Slip on the stirrup

2




Seam on the top of the foot

3




1st turn around the heel

4




2nd turn securing from the bottom

5



3rd turn securing from the top

6



Circular application

Accessories

Order in
France



FOAM BANDAGE Products for medical padding

Thickness	Width	Length	Reference
8 mm	8 cm	1 m	530809 1100
	10 cm		531009 1100
12 mm	10 cm	1 m	541009 1100



BREATHABLE

limits the risks of maceration thanks to the open cell foam⁽¹⁾



VARICO Products for medical padding

Optimisation of compression treatment by leveling out of anatomical relief.

Thanks to the open cell foam, the pad is permeable to air and helps preventing the risk of wound maceration⁽¹⁾.

Model	Dimensions	Packaging	Reference
0	63 mm x 87 mm	x 2 units	520009 1000
1	70 mm x 132 mm	x 1 unit	520009 1100
2	125 mm x 164 mm	x 1 unit	520009 1200
3	200 mm x 250 mm	x 1 unit	520009 1300



EASY TO USE

- can be washed and reused for the same patient
- can be cut to adapt to the size of the limb




TARGETED APPLICATION

- can be used in association to Thuasne® compressive bandages, Mobiderm® pads and bandages or Biflex Self Adjust
- allows the homogenization of compressive treatment by equalizing anatomical contours

Biflex Accessories

Order in
France



SELF-GRIPPING FASTENERS Ensure good support and prevent risk of injury

Packaging	Color	Reference
x 4 units	White	590309 1000



ELASTIC FASTENERS

Packaging	Color	Reference
x 4 units	Beige	590209 1000

(1) IFTH report n°180798R

Biflex Self Adjust

SINCE 1932

Adjustable compression system
for upper and lower limb

Indicated for the management of oedema (venous, lymphatic and caused by injury), chronic venous insufficiency and symptomatic treatment of lipoedema.

Biflex Self Adjust was developed to treat efficiently and with an optimal comfort all kinds of oedemas.

EFFICIENCY:

- The good pressure level on the limb and the reproductivity of the fitting for the HCP and the patient himself is ensured by the patented pressure adjuster directly integrated to the product.
- The design was chosen to have a complete covering of the limb and avoid oedemas migration.


COMFORT:

- The material was chosen to improve the comfort: the open cell foam allows moisture transfer.
- Easy donning and doffing thanks to the integrated loop and an optimised number of bandages.




Biflex Self Adjust

SINCE 1932




Ensure a minimum pressure level of 20 mmHg to 40 mmHg for lower limb and of 20 to 30 mmHg for upper limb and the fitting reproductivity*

Integrated pressure adjuster system on lower limb (calf and thigh).



Complete covering of the limb to avoid oedema migration

Overlapped bandages design and reversed closures.



Material allowing moisture transfer**

Open cell foam between the materials.



Easy donning and doffing

Integrated loop for calf, thigh and arm.

Minimum number of bandages.

How to use the pressure adjuster ?






The adjuster comes with 3 levels of pressure:

L low (light blue line), **M** medium (blue line) and **H** high (dark blue line).


Pressure is reached by pulling on the strap which slides the adjuster until the outside of the loop is over the desired level of pressure.

When putting it on and wearing it, ensure that no large folds appear.

If necessary, the segment pressure can be adjusted by adjusting the loops from bottom to top.



CALF




Size	cB	cC	cD
1	16 - 24 cm	24 - 35 cm	25 - 36 cm
	6.5 - 9.5"	9.5 - 14"	10 - 14"
2	20 - 28 cm	32 - 43 cm	32 - 43 cm
	8 - 11"	12.5 - 17"	12.5 - 17"
3	24 - 32 cm	40 - 51 cm	39 - 50 cm
	9.5 - 12.5"	15.5 - 20"	15.5 - 19.5"
4	28 - 36 cm	48 - 59 cm	46 - 57 cm
	11 - 14"	19 - 23"	18 - 22.5"

N

L


 $l < 36 \text{ cm} / 14"$ $l > 36 \text{ cm} / 14"$

KNEE



Size	cE
1	35 - 50 cm
	14 - 19.5"
2	50 - 65 cm
	19.5 - 25.5"
3	65 - 85 cm
	25.5 - 33.5"

THIGH




Size	cE	cG
1	34 - 45 cm	46 - 57 cm
	13.5 - 18"	18 - 22.5"
2	41 - 53 cm	53 - 65 cm
	16 - 21"	21 - 25.5"
3	49 - 60 cm	61 - 72 cm
	19.5 - 23.5"	24 - 28.5"
4	56 - 68 cm	68 - 80 cm
	22 - 27"	26.5 - 31.5"

N

L


 $l < 36 \text{ cm} / 14"$ $l > 36 \text{ cm} / 14"$

FOOT



Size	cE	cG
1	< 36 cm	19 - 24 cm
	< 14"	7.5 - 9.5"
2	< 36 cm	24 - 29 cm
	< 14"	9.5 - 11.5"
3	< 41 cm	29 - 34 cm
	< 16"	11.5 - 13.5"

ARM



Size	cG	cD	cC
1	25 - 32 cm	22 - 28 cm	15 - 19 cm
	10 - 12.5"	9 - 11"	6 - 7.5"
2	29 - 36 cm	24 - 30 cm	16 - 20 cm
	11.5 - 14"	9.5 - 12"	6.5 - 8"
3	33 - 40 cm	26 - 32 cm	17 - 21 cm
	13 - 15.5"	10 - 12.5"	6.5 - 8"
4	37 - 44 cm	28 - 34 cm	17 - 21 cm
	14.5 - 17.5"	11 - 13.5"	6.5 - 8"
5	41 - 48 cm	30 - 36 cm	18 - 22 cm
	16 - 19"	12 - 14"	7 - 8.5"
6	45 - 52 cm	31 - 37 cm	19 - 23 cm
	18 - 20.5"	12.5 - 14.5"	7.5 - 9"

HAND



Size	cE
1	14 - 19 cm
	5.5 - 7.5"
2	19 - 24 cm
	7.5 - 9.5"
3	24 - 29 cm
	9.5 - 11.5"

Sizes



DESIGNATION

Code article

Biflex Self Adjust calf

T1 regul	18300120100201
T1 tall	18300120100202
T2 regul	18300120200201
T2 tall	18300120200202
T3 regul	18300120300201
T3 tall	18300120300202
T4 regul	18300120400201
T4 tall	18300120400202



Biflex Self Adjust knee

T1	18200120100299
T2	18200120200299
T3	18200120300299



Biflex Self Adjust thigh

T1 regul	18250120100201
T1 tall	18250120100202
T2 regul	18250120200201
T2 tall	18250120200202
T3 regul	18250120300201
T3 tall	18250120300202
T4 regul	18250120400201
T4 tall	18250120400202



Biflex Self Adjust foot

T1	18100120100299
T2	18100120200299
T3	18100120300299



Biflex Self Adjust arm

T1	18350120100299
T2	18350120200299
T3	18350120300299
T4	18350120400299
T5	18350120500299
T6	18350120600299



Biflex Self Adjust hand

T1	18150120100299
T2	18150120200299
T3	18150120300299



Protective liner – socks

TU	31117 9011292 1
----	-----------------

*Thuasne® doesn't only offer compression,
discover our other fields of expertise.*



Orthopaedics

Ranges of belts and braces:
neck, shoulder, wrist, thumb,
elbow, knee, ankle, foot...



Lymphoedema and scars

Medical compression devices
for the treatment of lymphoedema
and severe burns.



Post breast cancer SILIMA®

External breast prostheses,
lingerie, swimwear.



Compression THERAFIRM®

Range of US standard
compressive garments.



Prosthetic and orthotic KNIT-RITE®

Range of textile solutions for prosthetics,
orthotics and other needs.

Ask for our catalogues

WINGS FOR YOUR HEALTH

Thuasne®

120, rue Marius Aufan – CS 10032
92309 Levallois-Perret Cedex – France
Tel. + 33 1 41 05 92 92 – Fax + 33 1 41 05 89 91
e-mail: sales.export@thuasne.fr

www.thuasne.com



Venoflex products: medical compression socks, stockings or tights for women and men.
Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.
All the medical devices mentioned on this document are CE marked according to the European Regulation 2017/745 on medical devices.
Products mentioned in this document are CE class I devices. Please contact Thuasne® should you need any additional information on devices classification.
Please read carefully the instructions for use, indications and contraindications of the products.

We reserve the right to modify our products as shown and described in this catalogue, without previous notice.
Non-contractual photographs.
Please check the availability of products according to the countries.
Ref.: 0117 092 100 (2026-02) - 2511391.
Photos: Studio Caterin - Daniel Reinhold, Hendrik Nennecke - StockLite/Shutterstock.com - Thuasne®.

