



Cicatrex

Please send your order to
export.m2m@thuasne.de

Company address

Contact person _____

Tel. number _____

Date _____

☐ Quotation by email:

☐ Order

Shipping address

Patient

Last name: _____ ☐ m

First name: _____ ☐ f

Date of birth: _____

Address: _____

☐ New supply

☐ Reorder, ref. No.:

☐ Reorder with changes

Internal notes

Quantity	Model	Material	Color	Seams with a multi-colored yarn	Iron-on patches
☐ Piece	☐ Special arm* ☐ Vest ☐ Jacket	☐ Nature knitted ☐ Nature woven ☐ Skin ☐ Rigid	☐ Beige ☐ Black ☐ Pink ☐ Mint ¹ (¹ not available in Skin, ² only available in Skin, ³ only available in Rigid/Nature knitted)	☐ Jeans ☐ Trend ² ☐ Dark blue ³	☐ Cat ☐ Owl ☐ Helicopter ☐ Spider ☐ Animals ☐ Flower girl ☐ Butterflies ☐ Flowers ☐ Bunny ☐ Football ☐ Cats ☐ Pirate ☐ Flag ☐ Unicorn

(* only available in Skin/Nature woven)

Special features (please enclose or draw in photos)			
☐ Amputation ☐ Pressure pad	☐ Memory foam ☐ Silon-TEX®	☐ Additional opening ☐ Soft lining fabric	☐ ARION donning aid size: _____

Thorax

Right side

- 13 _____ Right shoulder length
- 14 _____ Right underarm circumference
- 16 _____ Right shoulder cap
- 15 _____ Right upper arm circumference
- 18 _____ Length to mid-upper arm
- 17 _____ Mid-upper arm circumference (widest point)
- 20 _____ Length to elbow
- 19 _____ Elbow circumference
- 22 _____ Length to mid-forearm
- 21 _____ Circumference mid-forearm (widest point)
- 24 _____ Length to wrist
- 23 _____ Wrist circumference

- 1 _____ Neck circumference over 7th cervical vertebra and top edge of sternum
- 2 _____ Half underarm front circumference, inhaled
- 4 _____ Length from upper edge of sternum to underarm
- 3 _____ Circumferential underarm measurement, inhaled
- 6 _____ Length from upper edge of sternum to underbust
- 5 _____ Underbust circumference, inhaled
- 8 _____ Length from upper edge of sternum to waist
- 7 _____ Waist circumference
- 10 _____ Length from upper edge of sternum to end of the bandage
- 9 _____ Pelvic circumference
- 37 _____ Length 7th cervical vertebra to bandage end

Left side

- 12 _____ Left underarm circumference
- 11 _____ Left shoulder length
- 25 _____ Circumference upper arm left
- 26 _____ Left shoulder cap
- 27 _____ Mid-upper arm circumference (widest point)
- 28 _____ Length to mid-upper arm
- 29 _____ Elbow circumference
- 30 _____ Length to elbow
- 31 _____ Circumference mid-forearm (widest point)
- 32 _____ Length to mid-forearm
- 33 _____ Length to wrist
- 34 _____ Wrist circumference

Right arm

Armpit compression ☐ yes ☐ no

Fastener

- Length _____ cm
from: _____ to: _____
- ☐ Zipper
☐ sewn in ☐ sewn on
- ☐ Hooks and eyes
- ☐ Thumb side, front
- ☐ Little finger side, back
- ☐ Inside
- ☐ Outside
- ☐ Deficit of stretching (anatomical elbow): _____°
- ☐ Silicone dot band
☐ 3 cm ☐ 5 cm
- ☐ Elastic tape
☐ 3 cm ☐ 5 cm

- ☐ Neckline front
_____ cm from sternum
- ☐ Neckline back
_____ cm from 7 cervical whirl
- ☐ Collar
- Upper circumference: _____ cm
- Lower circumference: _____ cm
- Front height: _____ cm
- Rear height: _____ cm
- Sternum to lower neck _____ cm
- ☐ Cups
- ☐ with compression ☐ without compression
- right: A B C D E F
- left: A B C D E F
- ☐ Breast form pocket ☐ right ☐ left

Thorax

Fastener


- ☐ Zipper
- ☐ sewn in ☐ sewn on
- ☐ Hooks and eyes
- Position ☐ chest ☐ back
- ☐ Crotch strap removable _____ cm long
- ☐ Crotch strap fixed on one side
- ☐ Silicone dot band
☐ 3 cm ☐ 5 cm
- ☐ Elastic tape (without silicone)
☐ 3 cm ☐ 5 cm
- ☐ Waistband
- ☐ Velcro fastener to pants

Left arm

Armpit compression ☐ yes ☐ no

Fastener

- Length _____ cm
from: _____ to: _____
- ☐ Zipper
☐ sewn in ☐ sewn on
- ☐ Hooks and eyes
- ☐ Thumb side, front
- ☐ Little finger side, back
- ☐ Inside
- ☐ Outside
- ☐ Deficit of stretching (anatomical elbow): _____°
- ☐ Silicone dot band
☐ 3 cm ☐ 5 cm
- ☐ Elastic tape
☐ 3 cm ☐ 5 cm



Instructions for the accurate measurement of a Cicatrex jacket

- For precise measurement, use a flexible measuring tape and a body marker.
- All measurements should be taken while the patient is in a natural, relaxed standing posture.
- If possible, outline the jacket on the body sketch and shade in the scarred areas.
- Pressure pads, Silon-Tex® elements, additional openings, and any other measurements deemed relevant by the person taking the measurements must also be indicated in the sketch.
- Sleeve length measurements are to be taken starting from the acromion. Always measure along the outer side of the arm, passing over the lateral epicondyle.

Key to Symbols and Abbreviations

- ⌒ = Semi-circumferential measurement
- = Circumferential measurement
- = Length measurement
- Z = Zipper